



# JBL CristalProfi Pulling Aid Impell. Bearing&Brush

Pulling aid and cleaning brush for impeller bearings

Suitable for:   



- Practical auxiliary tool to maintain and to improve the performance (JBL CristalProfi e: 401, 402, 701, 702, 901, 902, 1501, 1502, 1901, 1902)
- Easy to use: remove the bearing bush at the end of the ceramic axis with the pulling aid. Clean the well with the brush
- Diameters: brush 17 mm, pulling aid 4.5 mm
- Cleans areas which are difficult to reach perfectly. The lower impeller bearing (rubber sleeve) in external filters is practically inaccessible without using the correct tool
- Package contents: pulling aid and cleaning brush for impeller bearings



# JBL CristalProfi Pulling Aid Impell. Bearing&Brush


## Product information

### Clean and healthy water

Plant and food remains and metabolic substances lead to a deterioration of the quality of the water in the aquarium. Good water quality is necessary for healthy fish and plants. This can be achieved with various filters. The filters suck in the aquarium water and remove the pollution and waste materials from the water. Filters also offer an ideal habitat for bacteria which decompose harmful substances.

### Maintenance of the pump performance

The pulling aid helps you pull the bearing bush out of the external filter at the end of the ceramic axis. Without a pulling aid it is difficult to reach. With the help of the brush the well, in which the impeller rotates, can be cleaned perfectly. Cleaning the impeller well improves the performance and restores the original performance of the pump.


Further information	
FAQ	✓
Blog	✓
Press	✓
Laboratory/calculator	✗
Worth reading	✓
Spare parts	✗
Video	✓
GarantiePlus	✗
Instructions	✓
QR code	



# JBL CristalProfi Pulling Aid Impell. Bearing&Brush



## Product details

Article data	
Product name	JBL CP pulling aid impell. bearing & brush
Art. No.	6011400
EAN number	4014162601148
EAN as barcode	
Content	1 kit
-	-
Expiry months	-
RRP incl. VAT	8,46 €
Base price	-
Nominal filling quantity	-
Base quantity	1
Gross weight	31 g
Net weight	11 g
Weight change	1000

Disposal	
Product name	JBL CP pulling aid impell. bearing & brush
Art. No.	6011400
Green dot	✓
Group electronic waste	-
Disposal weight	-
Battery type	-
Battery return	-
Battery rechargeable	-
Disposal weight battery	-
Non-returnable glass	-
PPK	10 g
Plastic small	10 g
Plastic large	-
Disposal weight metal	0 g

Features	
Product name	JBL CP pulling aid impell. bearing & brush
Art. No.	6011400





Features	
Animal species	Armored catfish, Arowana, Axolotl, Barbels, Bettas, Bichirs/ reedfish, Blowfish, Catfish, Cichlids (South America), Clawed frogs, Clawed shrimps, Corals, Crayfish, Crustaceans, Danions, Discus, Dwarf crayfish, Dwarf shrimps, Flowerhorn, Gill maggots, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, blood parrot cichlids, freshwater butterflyfish, turtle
Animal size	For all animal sizes
Animal age group	All aquarium fish
Volume of habitat	For all aquariums
Material	plastic (ABS), metal (steel)
Food type	-
Colour	black / white / silver
Dosage	-
Transport conditions	-



# JBL CristalProfi Pulling Aid Impell. Bearing&Brush

Electronic label / illuminant	
Product name	JBL CP pulling aid impell. bearing & brush
Art. No.	6011400
Ambient temperature	-
Start time	-
Mercury	-
Tube length	-
Service life	-
Lumen	-
CRI value	-
Dimmable	-
Switching cycles	-
PAR value	-
Energy efficiency class	-
UV-A	-
UV-B	-
UV-C	-
Colour temperature	-
Base designation	-

Technical data	
Product name	JBL CP pulling aid impell. bearing & brush
Art. No.	6011400
Range in litres	-
Range from - to	-
Range in days	-
Range tank length	-
Output in watts	-
Output per hour	-
Output per day	-
Height	-
Length	-
Width	-
Diameter	-
Voltage	-
For	-
T8 26mm (watt)	-
T5 16mm (watt)	-
Size	-
Content for	-
Filter container volume	-
Volume filter media	-
Hose connections pressure/out	-





Technical data	
Hose connections suction/in	-
Delivery head	-





## JBL CristalProfi Pulling Aid Impell. Bearing&Brush

Food type	-
Sub product type	-
Dosing	-





## Additional information for the specialist trade sector

Article data	
Product name	JBL CP pulling aid impell. bearing & brush
Art. No.	6011400
VAT	20%
Sales unit (SU)	3
Volume packaging	0.4l
Dimensions (l/w/h)	23 mm/96 mm/218 mm
Layer	360
Pallet	2160
Category of products	3
Customs tariff	96039099
Country of origin	CN
Type of packaging	Skin / blister card packaging

PU 1 data	
Product name	JBL CP pulling aid impell. bearing & brush
Art. No.	6011400
PU 1 material	cardboard gr
PU 1 weight	90 g
PU 1 length	100 mm
PU 1 width/depth	70 mm
PU 1 height	220 mm

PU 2 data	
Product name	JBL CP pulling aid impell. bearing & brush
Art. No.	6011400
PU 2 material	cardboard gr
PU 2 weight	320 g
PU 2 length	400 mm
PU 2 width/depth	300 mm
PU 2 height	300 mm

Trade data	
Product name	JBL CP pulling aid impell. bearing & brush
Art. No.	6011400
Till receipt text	CP e Zughilfe
Shelf placement	-

