



JBL CristalProfi Pulling Aid Impell. Bearing&Brush

Pulling aid and cleaning brush for impeller bearings

Suitable for:   



- Practical auxiliary tool to maintain and to improve the performance (JBL CristalProfi e: 401, 402, 701, 702, 901, 902, 1501, 1502, 1901, 1902)
- Easy to use: remove the bearing bush at the end of the ceramic axis with the pulling aid. Clean the well with the brush
- Diameters: brush 17 mm, pulling aid 4.5 mm
- Cleans areas which are difficult to reach perfectly. The lower impeller bearing (rubber sleeve) in external filters is practically inaccessible without using the correct tool
- Package contents: pulling aid and cleaning brush for impeller bearings



JBL CristalProfi Pulling Aid Impell. Bearing&Brush




Product information

Clean and healthy water

Plant and food remains and metabolic substances lead to a deterioration of the quality of the water in the aquarium. Good water quality is necessary for healthy fish and plants. This can be achieved with various filters. The filters suck in the aquarium water and remove the pollution and waste materials from the water. Filters also offer an ideal habitat for bacteria which decompose harmful substances.

Maintenance of the pump performance

The pulling aid helps you pull the bearing bush out of the external filter at the end of the ceramic axis. Without a pulling aid it is difficult to reach. With the help of the brush the well, in which the impeller rotates, can be cleaned perfectly. Cleaning the impeller well improves the performance and restores the original performance of the pump.

| Further information | |
|-----------------------|---|
| FAQ | ✓ |
| Blog | ✓ |
| Press | ✓ |
| Laboratory/calculator | ✗ |
| Worth reading | ✓ |
| Spare parts | ✗ |
| Video | ✓ |
| GarantiePlus | ✗ |
| Instructions | ✓ |
| QR code |  |






JBL CristalProfi Pulling Aid Impell. Bearing&Brush



Product details

| Article data | |
|--------------------------|---|
| Product name | JBL CP pulling aid impell. bearing & brush |
| Art. No. | 6011400 |
| EAN number | 4014162601148 |
| EAN as barcode |  |
| Content | 1 kit |
| - | - |
| Expiry months | - |
| RRP incl. VAT | 8,39 € |
| Base price | - |
| Nominal filling quantity | - |
| Base quantity | 1 |
| Gross weight | 31 g |
| Net weight | 11 g |
| Weight change | 1000 |

| Disposal | |
|-------------------------|--|
| Product name | JBL CP pulling aid impell. bearing & brush |
| Art. No. | 6011400 |
| Green dot | ✓ |
| Group electronic waste | - |
| Disposal weight | - |
| Battery type | - |
| Battery return | - |
| Battery rechargeable | - |
| Disposal weight battery | - |
| Non-returnable glass | - |
| PPK | 10 g |
| Plastic small | 10 g |
| Plastic large | - |
| Disposal weight metal | 0 g |

| Features | |
|--------------|--|
| Product name | JBL CP pulling aid impell. bearing & brush |
| Art. No. | 6011400 |





| Features | |
|----------------------|--|
| Animal species | Armored catfish, Arowana, Axolotl, Barbels, Bettas, Bichirs/ reedfish, Blowfish, Catfish, Cichlids (South America), Clawed frogs, Clawed shrimps, Corals, Crayfish, Crustaceans, Danions, Discus, Dwarf crayfish, Dwarf shrimps, Flowerhorn, Gill maggots, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, blood parrot cichlids, freshwater butterflyfish, turtle |
| Animal size | For all animal sizes |
| Animal age group | All aquarium fish |
| Volume of habitat | For all aquariums |
| Material | plastic (ABS), metal (steel) |
| Food type | - |
| Colour | black / white / silver |
| Dosage | - |
| Transport conditions | - |



JBL CristalProfi Pulling Aid Impell. Bearing&Brush

| Electronic label / illuminant | |
|-------------------------------|--|
| Product name | JBL CP pulling aid impell. bearing & brush |
| Art. No. | 6011400 |
| Ambient temperature | - |
| Start time | - |
| Mercury | - |
| Tube length | - |
| Service life | - |
| Lumen | - |
| CRI value | - |
| Dimmable | - |
| Switching cycles | - |
| PAR value | - |
| Energy efficiency class | - |
| UV-A | - |
| UV-B | - |
| UV-C | - |
| Colour temperature | - |
| Base designation | - |

| Technical data | |
|-------------------------------|--|
| Product name | JBL CP pulling aid impell. bearing & brush |
| Art. No. | 6011400 |
| Range in litres | - |
| Range from - to | - |
| Range in days | - |
| Range tank length | - |
| Output in watts | - |
| Output per hour | - |
| Output per day | - |
| Height | - |
| Length | - |
| Width | - |
| Diameter | - |
| Voltage | - |
| For | - |
| T8 26mm (watt) | - |
| T5 16mm (watt) | - |
| Size | - |
| Content for | - |
| Filter container volume | - |
| Volume filter media | - |
| Hose connections pressure/out | - |





| Technical data | |
|-----------------------------|---|
| Hose connections suction/in | - |
| Delivery head | - |



JBL CristalProfi Pulling Aid Impell. Bearing&Brush

| | |
|------------------|---|
| Food type | - |
| Sub product type | - |
| Dosing | - |





Additional information for the specialist trade sector

| Article data | |
|----------------------|--|
| Product name | JBL CP pulling aid impell. bearing & brush |
| Art. No. | 6011400 |
| VAT | 19% |
| Sales unit (SU) | 3 |
| Volume packaging | 0.4l |
| Dimensions (l/w/h) | 23 mm/96 mm/218 mm |
| Layer | 360 |
| Pallet | 2160 |
| Category of products | 3 |
| Customs tariff | 96039099 |
| Country of origin | CN |
| Type of packaging | Skin / blister card packaging |

| PU 1 data | |
|------------------|--|
| Product name | JBL CP pulling aid impell. bearing & brush |
| Art. No. | 6011400 |
| PU 1 material | cardboard gr |
| PU 1 weight | 90 g |
| PU 1 length | 100 mm |
| PU 1 width/depth | 70 mm |
| PU 1 height | 220 mm |

| PU 2 data | |
|------------------|--|
| Product name | JBL CP pulling aid impell. bearing & brush |
| Art. No. | 6011400 |
| PU 2 material | cardboard gr |
| PU 2 weight | 320 g |
| PU 2 length | 400 mm |
| PU 2 width/depth | 300 mm |
| PU 2 height | 300 mm |

| Trade data | |
|-------------------|--|
| Product name | JBL CP pulling aid impell. bearing & brush |
| Art. No. | 6011400 |
| Till receipt text | CP e Zughilfe |
| Shelf placement | - |