



JBL PROSCAPE CLEANING GLOVE

Aquarium cleaning glove

Suitable for:



- Effortless removal of algae growth and deposits: aquarium glove with metal threads for the cleaning of glass panes and objects
- Easy handling: put on glove, clean by rubbing on the desired area
- Cleans all corners and angles, removes algae from plant leaves and decoration, cleans technical items and pumps
- Care of glove: rinse in lukewarm water and dry out of direct sunlight
- Package contents: 1x aquarium glove



You may also be interested in

You can find a complete overview here: <https://www.jbl.de/qr/61379>



JBL Aqua-T Handy
Glass pane cleaner with stainless steel blade



JBL Aqua-T Triumph
Glass pane cleaner with spare blades and rubber scraper



JBL FLOATY acrylic/glass
Floating acrylic pane cleaning magnet



JBL Algae Magnet
Magnetic glass cleaner for aquarium panes





JBL PROSCAPE CLEANING GLOVE



Product information

Maintenance and cleaning


Food and plant remains and metabolic products form impurities in an aquarium. Cleaning your tank and accessories and carrying out regular partial water changes is good for your aquarium dwellers and makes having an aquarium more fun for you.

Effortless cleaning

The JBL ProScape cleaning glove is the ideal helper for the fast and effortless removal of algae growth and deposits on aquarium panes, devices and decorations.

Care instructions

Rinse in lukewarm water and dry out of direct sun light.

| Further information | |
|-----------------------|---|
| FAQ | ✓ |
| Blog | ✓ |
| Press | ✓ |
| Laboratory/calculator | ✗ |
| Worth reading | ✓ |
| Spare parts | ✗ |
| Video | ✓ |
| GarantiePlus | ✗ |
| Instructions | ✓ |
| QR code |  |





JBL PROSCAPE CLEANING GLOVE

Product details

| Article data | | | |
|--------------------------|-----------------------------|---------------------------------|----------------------------------|
| Product name | JBL PROSCAPE CLEANING GLOVE | JBL PROSCAPE CLEANING GLOVE S/M | JBL PROSCAPE CLEANING GLOVE L/XL |
| Art. No. | 6137900 | 6138600 | 6138700 |
| EAN number | 4014162613790 | 4014162613868 | 4014162613875 |
| EAN as barcode | | | |
| Content | 1 piece(s) | - | - |
| - | - | - | - |
| Expiry months | - | - | - |
| RRP incl. VAT | 13,07 € | 13,74 € | 13,74 € |
| Base price | - | - | - |
| Nominal filling quantity | - | - | - |
| Base quantity | 1 | - | - |
| Gross weight | 30 g | 30 g | 30 g |
| Net weight | 26 g | 26 g | 26 g |
| Weight change | 1000 | - | - |

| Disposal | | | |
|-------------------------|-----------------------------|---------------------------------|----------------------------------|
| Product name | JBL PROSCAPE CLEANING GLOVE | JBL PROSCAPE CLEANING GLOVE S/M | JBL PROSCAPE CLEANING GLOVE L/XL |
| Art. No. | 6137900 | 6138600 | 6138700 |
| Green dot | ✓ | - | - |
| Group electronic waste | - | - | - |
| Disposal weight | - | - | - |
| Battery type | - | - | - |
| Battery return | - | - | - |
| Battery rechargeable | - | - | - |
| Disposal weight battery | - | - | - |
| Non-returnable glass | - | - | - |
| PPK | - | - | - |
| Plastic small | 4 g | - | - |
| Plastic large | - | - | - |
| Disposal weight metal | 0 g | 0 g | 0 g |

| Features | | | |
|--------------|-----------------------------|---------------------------------|----------------------------------|
| Product name | JBL PROSCAPE CLEANING GLOVE | JBL PROSCAPE CLEANING GLOVE S/M | JBL PROSCAPE CLEANING GLOVE L/XL |
| Art. No. | 6137900 | 6138600 | 6138700 |





| Features | | | |
|----------------------|--|--|--|
| Animal species | Armored catfish, Arowana, Axolotl, Barbels, Bettas, Bichirs/ reedfish, Blowfish, Catfish, Cichlids (South America), Clawed frogs, Clawed shrimps, Crayfish, Crustaceans, Danions, Discus, Dwarf crayfish, Dwarf shrimps, Flowerhorn, Gill maggots, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, blood parrot cichlids, freshwater butterflyfish, turtle | Armored catfish, Arowana, Axolotl, Barbels, Bettas, Bichirs/ reedfish, Blowfish, Catfish, Cichlids (South America), Clawed frogs, Clawed shrimps, Crayfish, Crustaceans, Danions, Discus, Dwarf crayfish, Dwarf shrimps, Flowerhorn, Gill maggots, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, blood parrot cichlids, freshwater butterflyfish, turtle | Armored catfish, Arowana, Axolotl, Barbels, Bettas, Bichirs/ reedfish, Blowfish, Catfish, Cichlids (South America), Clawed frogs, Clawed shrimps, Crayfish, Crustaceans, Danions, Discus, Dwarf crayfish, Dwarf shrimps, Flowerhorn, Gill maggots, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, blood parrot cichlids, freshwater butterflyfish, turtle |
| Animal size | For all animal sizes | For all animal sizes | For all animal sizes |
| Animal age group | All aquarium fish | All aquarium fish | All aquarium fish |
| Volume of habitat | For all aquariums | - | - |
| Material | acetate with lurex | acetate with lurex | acetate with lurex |
| Food type | - | - | - |
| Colour | silver / white | silver / white | silver / white |
| Dosage | - | - | - |
| Transport conditions | - | - | - |



JBL PROSCAPE CLEANING GLOVE

| Electronic label / illuminant | | | |
|-------------------------------|-----------------------------|---------------------------------|----------------------------------|
| Product name | JBL PROSCAPE CLEANING GLOVE | JBL PROSCAPE CLEANING GLOVE S/M | JBL PROSCAPE CLEANING GLOVE L/XL |
| Art. No. | 6137900 | 6138600 | 6138700 |
| Ambient temperature | - | - | - |
| Start time | - | - | - |
| Mercury | - | - | - |
| Tube length | - | - | - |
| Service life | - | - | - |
| Lumen | - | - | - |
| CRI value | - | - | - |
| Dimmable | - | - | - |
| Switching cycles | - | - | - |
| PAR value | - | - | - |
| Energy efficiency class | - | - | - |
| UV-A | - | - | - |
| UV-B | - | - | - |
| UV-C | - | - | - |
| Colour temperature | - | - | - |
| Base designation | - | - | - |

| Technical data | | | |
|-------------------------------|-----------------------------|---------------------------------|----------------------------------|
| Product name | JBL PROSCAPE CLEANING GLOVE | JBL PROSCAPE CLEANING GLOVE S/M | JBL PROSCAPE CLEANING GLOVE L/XL |
| Art. No. | 6137900 | 6138600 | 6138700 |
| Range in litres | - | - | - |
| Range from - to | - | - | - |
| Range in days | - | - | - |
| Range tank length | - | - | - |
| Output in watts | - | - | - |
| Output per hour | - | - | - |
| Output per day | - | - | - |
| Height | - | - | - |
| Length | - | - | - |
| Width | - | - | - |
| Diameter | - | - | - |
| Voltage | - | - | - |
| For | - | - | - |
| T8 26mm (watt) | - | - | - |
| T5 16mm (watt) | - | - | - |
| Size | - | - | - |
| Content for | - | - | - |
| Filter container volume | - | - | - |
| Volume filter media | - | - | - |
| Hose connections pressure/out | - | - | - |
| Hose connections suction/in | - | - | - |





| Technical data | | | |
|----------------|---|---|---|
| Delivery head | - | - | - |





JBL PROSCAPE CLEANING GLOVE

| | |
|------------------|---|
| Food type | - |
| Sub product type | - |
| Dosing | - |





Additional information for the specialist trade sector

| Article data | | | |
|----------------------|-----------------------------|---------------------------------|----------------------------------|
| Product name | JBL PROSCAPE CLEANING GLOVE | JBL PROSCAPE CLEANING GLOVE S/M | JBL PROSCAPE CLEANING GLOVE L/XL |
| Art. No. | 6137900 | 6138600 | 6138700 |
| VAT | 19% | 19% | 19% |
| Sales unit (SU) | 6 | 6 | 6 |
| Volume packaging | 0.36l | 0.36l | 0.36l |
| Dimensions (l/w/h) | 6 mm/150 mm/300 mm | 6 mm/150 mm/300 mm | 6 mm/150 mm/300 mm |
| Layer | 1152 | 1152 | 1152 |
| Pallet | 2304 | 2304 | 2304 |
| Category of products | 2 | 2 | 2 |
| Customs tariff | 73231000 | 73231000 | 73231000 |
| Country of origin | JP | JP | JP |
| Type of packaging | piece(s) | | |

| PU 1 data | | | |
|------------------|-----------------------------|---------------------------------|----------------------------------|
| Product name | JBL PROSCAPE CLEANING GLOVE | JBL PROSCAPE CLEANING GLOVE S/M | JBL PROSCAPE CLEANING GLOVE L/XL |
| Art. No. | 6137900 | 6138600 | 6138700 |
| PU 1 material | film gr | - | - |
| PU 1 weight | 15 g | - | - |
| PU 1 length | 150 mm | - | - |
| PU 1 width/depth | 36 mm | - | - |
| PU 1 height | 300 mm | - | - |

| PU 2 data | | | |
|------------------|-----------------------------|---------------------------------|----------------------------------|
| Product name | JBL PROSCAPE CLEANING GLOVE | JBL PROSCAPE CLEANING GLOVE S/M | JBL PROSCAPE CLEANING GLOVE L/XL |
| Art. No. | 6137900 | 6138600 | 6138700 |
| PU 2 material | cardboard gr | - | - |
| PU 2 weight | 365 g | - | - |
| PU 2 length | 385 mm | - | - |
| PU 2 width/depth | 250 mm | - | - |
| PU 2 height | 270 mm | - | - |

| Trade data | | | |
|-------------------|-----------------------------|---------------------------------|----------------------------------|
| Product name | JBL PROSCAPE CLEANING GLOVE | JBL PROSCAPE CLEANING GLOVE S/M | JBL PROSCAPE CLEANING GLOVE L/XL |
| Art. No. | 6137900 | 6138600 | 6138700 |
| Till receipt text | Cleaning Glove | - | - |
| Shelf placement | - | - | - |