




# JBL CRISTALPROFI UNIBLOC e

Bio-filter foam insert for CristalProfi e

Suitable for:



- Clean and healthy water: bio-filter foam insert for the middle baskets of CristalProfi e
- Use: insert into the middle filter basket
- Remove round middle section for special filter materials
- Pore size especially suited for bio-filtering
- Package contents: 1 bio-filter foam 25 ppi insert in pack of 2

You may also be interested in

You can find a complete overview here: <https://www.jbl.de/qr/60161>



**JBL Symec micro**  
 Microfibre fleece against water cloudiness



**JBL Symec Filter floss**  
 Filter wool to remove all types of cloudiness



**JBL Symec XL**  
 Coarse filter wool to remove all types of water cloudiness



**JBL Symec VL**  
 Filter wool fleece to remove water cloudiness





# JBL CRISTALPROFI UNIBLOC e



## Product information

### Quality of foam:

Foams used in aquaristics with live animals must always be subjected to a pollutant test, as most plastics release problematic plasticisers into the water. All the plastics used by JBL have passed this test and do NOT release any harmful substances into the water.

### Cleaning of plastics:


For many months it is sufficient to rinse the dirty filter sponges under lukewarm running water and then inoculate them with filter bacteria (JBL FilterStart) before using them again. However, as mechanical cleaning NEVER leads to complete cleaning and dirt will remain in the filter sponge, it is recommended you replace the filter sponges about every 6 months.

If you place a freshly cleaned filter sponge in a glass of clean water overnight, you will be able to smell how much dirt is still present in the supposedly clean sponge the next morning!

### Inoculating filter material with cleansing bacteria:

Initial colonisation of filter material with cleansing bacteria takes about 1-2 weeks in freshwater and twice as long in seawater. If the filter material was already in use, the time for a "recolonisation" is doubled, as recolonisation does not take place so easily on dead bacterial lawns. The addition of concentrated cleaning bacteria such as in JBL FilterStart shortens the biological start-up time to about 24 hours.

Another reason is to fill the existing settlement space in the filter material with beneficial bacteria so that pathogenic germs find no place to spread. This has now been scientifically proven.

Further information	
FAQ	✓
Blog	✓
Press	✗
Laboratory/calculator	✗
Worth reading	✓
Spare parts	✗
Video	✓
GarantiePlus	✗
Instructions	✓
QR code	





# JBL CRISTALPROFI UNIBLOC e



## Product details

Article data		
Product name	JBL CRISTALPROFI UNIBLOC e4/7/901-2	JBL CRISTALPROFI UNIBLOC e15/1901-2
Art. No.	6016100	6016200
EAN number	4014162601612	4014162601629
EAN as barcode		
Content	2 sets	2 sets
For CPe	4-7-901,2	15-1901,2
Expiry months	-	-
RRP incl. VAT	11,77 €	17,22 €
Base price	-	-
Nominal filling quantity	-	-
Base quantity	1	1
Gross weight	100 g	126 g
Net weight	24 g	36 g
Weight change	1000	1000

Disposal		
Product name	JBL CRISTALPROFI UNIBLOC e4/7/901-2	JBL CRISTALPROFI UNIBLOC e15/1901-2
Art. No.	6016100	6016200
Green dot	✓	✓
Group electronic waste	-	-
Disposal weight	-	-
Battery type	-	-
Battery return	-	-
Battery rechargeable	-	-
Disposal weight battery	-	-
Non-returnable glass	-	-
PPK	76 g	90 g
Plastic small	-	-
Plastic large	-	-
Disposal weight metal	0 g	0 g

Features		
Product name	JBL CRISTALPROFI UNIBLOC e4/7/901-2	JBL CRISTALPROFI UNIBLOC e15/1901-2
Art. No.	6016100	6016200





Features		
Animal species	Armored catfish, Arowana, Axolotl, Barbels, Bettas, Bichirs/ reedfish, Blowfish, Catfish, Cichlids (South America), Clawed frogs, Clawed shrimps, Corals, Crayfish, Crustaceans, Danions, Discus, Dwarf crayfish, Dwarf shrimps, Flowerhorn, Gill maggots, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, blood parrot cichlids, freshwater butterflyfish, turtle	Armored catfish, Arowana, Axolotl, Barbels, Bettas, Bichirs/ reedfish, Blowfish, Catfish, Cichlids (South America), Clawed frogs, Clawed shrimps, Corals, Crayfish, Crustaceans, Danions, Discus, Dwarf crayfish, Dwarf shrimps, Flowerhorn, Gill maggots, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, blood parrot cichlids, freshwater butterflyfish, turtle
Animal size	For all animal sizes	For all animal sizes
Animal age group	All aquarium fish	All aquarium fish
Volume of habitat	> 300 L	800 L
Material	polyether	polyether
Food type	-	-
Colour	blue	blue
Dosage	-	-
Transport conditions	-	-





# JBL CRISTALPROFI UNIBLOC e

Electronic label / illuminant		
Product name	JBL CRISTALPROFI UNIBLOC e4/7/901-2	JBL CRISTALPROFI UNIBLOC e15/1901-2
Art. No.	6016100	6016200
Ambient temperature	-	-
Start time	-	-
Mercury	-	-
Tube length	-	-
Service life	-	-
Lumen	-	-
CRI value	-	-
Dimmable	-	-
Switching cycles	-	-
PAR value	-	-
Energy efficiency class	-	-
UV-A	-	-
UV-B	-	-
UV-C	-	-
Colour temperature	-	-
Base designation	-	-

Technical data		
Product name	JBL CRISTALPROFI UNIBLOC e4/7/901-2	JBL CRISTALPROFI UNIBLOC e15/1901-2
Art. No.	6016100	6016200
Range in litres	-	-
Range from - to	-	-
Range in days	-	-
Range tank length	-	-
Output in watts	-	-
Output per hour	-	-
Output per day	-	-
Height	-	-
Length	-	-
Width	-	-
Diameter	-	-
Voltage	-	-
For	-	-
T8 26mm (watt)	-	-
T5 16mm (watt)	-	-
Size	-	-
Content for	-	-
Filter container volume	-	-
Volume filter media	-	-
Hose connections pressure/out	-	-
Hose connections suction/in	-	-





Technical data		
Delivery head	-	-





# JBL CRISTALPROFI UNIBLOC e

Food type	-
Sub product type	-
Dosing	-





## Additional information for the specialist trade sector

Article data		
Product name	JBL CRISTALPROFI UNIBLOC e4/7/901-2	JBL CRISTALPROFI UNIBLOC e15/1901-2
Art. No.	6016100	6016200
VAT	19%	19%
Sales unit (SU)	1	1
Volume packaging	2.03l	2.6l
Dimensions (l/w/h)	70 mm/165 mm/176 mm	75 mm/175 mm/197 mm
Layer	120	96
Pallet	360	384
Category of products	2	2
Customs tariff	39269097	39269097
Country of origin	CN	CN
Type of packaging	Faltschachtel/Karton	Faltschachtel/Karton

PU 1 data		
Product name	JBL CRISTALPROFI UNIBLOC e4/7/901-2	JBL CRISTALPROFI UNIBLOC e15/1901-2
Art. No.	6016100	6016200
PU 1 material	film gr	film gr
PU 1 weight	3.7 g	3.8 g
PU 1 length	480 mm	520 mm
PU 1 width/depth	70 mm	75 mm
PU 1 height	180 mm	205 mm

PU 2 data		
Product name	JBL CRISTALPROFI UNIBLOC e4/7/901-2	JBL CRISTALPROFI UNIBLOC e15/1901-2
Art. No.	6016100	6016200
PU 2 material	cardboard gr	cardboard gr
PU 2 weight	650 g	720 g
PU 2 length	500 mm	542 mm
PU 2 width/depth	368 mm	335 mm
PU 2 height	362 mm	430 mm

Trade data		
Product name	JBL CRISTALPROFI UNIBLOC e4/7/901-2	JBL CRISTALPROFI UNIBLOC e15/1901-2
Art. No.	6016100	6016200
Till receipt text	UB e700/e900	UB e1500
Shelf placement	-	-

