Test case with most important water tests + iron test

Suitable for: < 🖈 🏏 🐕 🐋

CARE

PFLEGE **ENTRETIEN** 







- Monitors the water values in aquariums simply and reliably. Determines: pH value, carbonate hardness, iron, nitrite, nitrate, CO2 (calculation via table)
- Easy to use: syringe for exact measurement of the water quantity, detailed instructions for use
- Laboratory comparator system to compensate the inherent water colouring, analysis record sheet, instructions for problem solving, CO2 table
- Child-proof reagent bottles, water-proof plastic case, refillable kit
- Package contents: 1 test case for planted aquariums, 5 reagents, 3 glass vials, syringe, dosing spoon, comparator block, plastic cuvette, colour charts, analysis record sheets, instructions



You may also be interested in You can find a complete overview here: https://www.jbl.de/gr/25500



JBL pH Test 3.0-10.0 Quick test to determine the acidity



JBL pH Test 6.0-7.6 pH test for freshwater aquariums in the range 6.0-7.6



JBL pH Test 7.4-9.0 pH test in aquariums and ponds in the range 7.4-9.0



JBL GH Test Quick test to determine the total hardness





Accessories



JBL Aquarium Thermometer Float Aquarium thermometer











JBL component set for water tests



JBL Testlab empty + inserts

JBL assortment Test Colour Charts Colour charts for JBL water tests







#### Product information

Healthy aquariums with conditions close to nature The right water values are dependent on the fish stock and the plants in the aquarium. Even if the water looks clear it can be contaminated. With bad values diseases or algae can appear in the aquarium. To maintain a healthy aquarium with conditions close to nature it is important to check and adapt the water values regularly. With the JBL Test Combi Set you can determine the 5 most important water values in your aquarium.

Precise measurement of the following values: pH test: acidity of the water from 3.0 to 10 KH test: pH stability of the water (carbonate hardness) Fe test: iron value to control the fertilisation NO2 test: toxic nitrogen compound NO3 test: reason for algae growth

To determine the CO2 content you just need to measure the pH and KH values. From both values arises the CO2 content in the water, which then can be referenced from the enclosed table.

You can analyse your measured water values in the JBL Online Laboratory free of charge and receive recommendations.

| Further information   |   |
|-----------------------|---|
| FAQ                   | ~ |
| Blog                  | ✓ |
| Press                 | ~ |
| Laboratory/calculator | × |
| Worth reading         | ~ |
| Spare parts           | ✓ |
| Video                 | × |
| GarantiePlus          | × |
| Instructions          | ~ |
| QR code               |   |





| Article data             |                            |
|--------------------------|----------------------------|
| Product name             | JBL Test Combi Set Plus Fe |
| Art. No.                 | 2550000                    |
| EAN number               | 4014162255006              |
| EAN as barcode           |                            |
| Content                  | 6 test kits                |
| -                        | -                          |
| Expiry months            | 24                         |
| RRP incl. VAT            | -                          |
| Base price               | -                          |
| Nominal filling quantity | -                          |
| Base quantity            | 1                          |
| Gross weight             | 524 g                      |
| Net weight               | 524 g                      |
| Weight change            | 1000                       |

| Disposal                |                            |
|-------------------------|----------------------------|
| Product name            | JBL Test Combi Set Plus Fe |
| Art. No.                | 2550000                    |
| Green dot               | <ul> <li>✓</li> </ul>      |
| Group electronic waste  | -                          |
| Disposal weight         | -                          |
| Battery type            | -                          |
| Battery return          | -                          |
| Battery rechargeable    | -                          |
| Disposal weight battery | -                          |
| Non-returnable glass    | 13 g                       |
| РРК                     | -                          |
| Plastic small           | 62 g                       |
| Plastic large           | -                          |
| Disposal weight metal   | 0 g                        |

| Features       |  |
|----------------|--|
| Product name   | JBL Test Combi Set Plus Fe   |
| Art. No.       | 2550000  |
| Animal species | Arowana, Axolotl, Barbels, Bettas,<br>Bichirs/reedfish, Blowfish, Catfish,<br>Cichlids (South America), Corals,<br>Crayfish, Danions, Discus, Dwarf<br>shrimps, Flowerhorn, Gobies,<br>Goldfish, Gouramis, Guppy,<br>Juvenile fish, Killifish, Livebearers,<br>Loaches, Mussels, Panchaxes,<br>Rainbowfish, Snails, Spiny eels,<br>Terrapins, Tetra, Tropical<br>terrapins, Veiltails, freshwater<br>butterflyfish, turtle |







| Features             |                        |
|----------------------|------------------------|
| Animal size          | For all animal sizes   |
| Animal age group     | All aquarium fish      |
| Volume of habitat    | For all aquariums      |
| Material             | -                      |
| Food type            | -                      |
| Colour               | multicoloured          |
| Dosage               | see instruction manual |
| Transport conditions | -                      |





| Electronic label / illuminant |                            |
|-------------------------------|----------------------------|
| Product name                  | JBL Test Combi Set Plus Fe |
| Art. No.                      | 2550000                    |
| Ambient temperature           | -                          |
| Start time                    | -                          |
| Mercury                       | -                          |
| Tube length                   | -                          |
| Service life                  | -                          |
| Lumen                         | -                          |
| CRI value                     | -                          |
| Dimmable                      | -                          |
| Switching cycles              | -                          |
| PAR value                     | -                          |
| Energy efficiency class       | -                          |
| UV-A                          | -                          |
| UV-B                          | -                          |
| UV-C                          | -                          |
| Colour temperature            | -                          |
| Base designation              | -                          |
|                               |                            |

| Technical data                |                            |
|-------------------------------|----------------------------|
| Product name                  | JBL Test Combi Set Plus Fe |
| Art. No.                      | 2550000                    |
| Range in litres               | -                          |
| Range from - to               | -                          |
| Range in days                 | -                          |
| Range tank length             | -                          |
| Output in watts               | -                          |
| Output per hour               | -                          |
| Output per day                | -                          |
| Height                        | -                          |
| Length                        | -                          |
| Width                         | -                          |
| Diameter                      | -                          |
| Voltage                       | -                          |
| For                           | -                          |
| T8 26mm (watt)                | -                          |
| T5 16mm (watt)                | -                          |
| Size                          | -                          |
| Content for                   | -                          |
| Filter container volume       | -                          |
| Volume filter media           | -                          |
| Hose connections pressure/out | -                          |
| Hose connections suction/in   | -                          |
| Delivery head                 | -                          |





| Food type        | -                      |
|------------------|------------------------|
| Sub product type | -                      |
| Dosing           | see instruction manual |



Safety information

| - safety information in accordance with GHS |  |  |
|---|--|--|
| Signal word                                 | Danger   |  |
| Hazard symbol                               |  |  |
| Risk phrase                                 | <ul> <li>H225: Highly flammable liquid and vapour.</li> <li>H412: Harmful to aquatic life with long lasting effects.</li> <li>H335: May cause respiratory irritation.</li> <li>H341: Suspected of causing genetic defects.</li> <li>H317: May cause an allergic skin reaction.</li> <li>H314: Causes severe skin burns and eye damage.</li> <li>H302: Harmful if swallowed.</li> <li>H261: In contact with water releases flammable gases.</li> </ul>  |  |
| safety note                                 | <ul> <li>P101: If medical advice is needed, have product container or label at hand.</li> <li>P102: Keep out of reach of children.</li> <li>P103: Read label before use.</li> <li>P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.</li> <li>P301+P310a: IF SWALLOWED: Immediately call a POISON CENTER/ doctor.</li> <li>P303+P361+P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.</li> <li>P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.</li> <li>P405: Store locked up.</li> <li>P501a: Dispose of contents/container in accordance with local/regional/ national/international regulations.</li> </ul> |  |

| Further safety information |              |
|----------------------------|--------------|
| Safety instructions        | $\checkmark$ |
| QR code                    |              |







#### Additional information for the specialist trade sector

| Article data         |                            |
|----------------------|----------------------------|
| Product name         | JBL Test Combi Set Plus Fe |
| Art. No.             | 2550000                    |
| VAT                  | 19%                        |
| Sales unit (SU)      | 1                          |
| Volume packaging     | 3.13I                      |
| Dimensions (l/w/h)   | 52 mm/266 mm/226 mm        |
| Layer                | 120                        |
| Pallet               | 360                        |
| Category of products | 1                          |
| Customs tariff       | 38220000                   |
| Country of origin    | DE                         |
| Type of packaging    | Faltschachtel/Karton       |

| PU 1 data        |                            |
|------------------|----------------------------|
| Product name     | JBL Test Combi Set Plus Fe |
| Art. No.         | 2550000                    |
| PU 1 material    | cardboard gr               |
| PU 1 weight      | 500 g                      |
| PU 1 lengh       | 540 mm                     |
| PU 1 width/depth | 260 mm                     |
| PU 1 height      | 470 mm                     |

| PU 2 data        |                            |
|------------------|----------------------------|
| Product name     | JBL Test Combi Set Plus Fe |
| Art. No.         | 2550000                    |
| PU 2 material    | None                       |
| PU 2 weight      | -                          |
| PU 2 lengh       | -                          |
| PU 2 width/depth | -                          |
| PU 2 height      | -                          |

| Trade data        |                            |
|-------------------|----------------------------|
| Product name      | JBL Test Combi Set Plus Fe |
| Art. No.          | 2550000                    |
| Till receipt text | Test Combi Set             |
| Shelf placement   | -                          |

