



Filter wool to remove all types of cloudiness









- Fast removal of all kinds of cloudiness in freshwater and marine aquariums: quality filter wool for aquarium filters
- Easy to use: select quantity of wool (suitable for all filters) as required. Insert as last filter stage. Removal required when water flow rate decreases.
- Crystal-clear water: filters off all types of cloudiness and particles. Although rinsing is possible we recommend exchanging it instead
- No release of pollutants and fibres fully synthetic, residue-free filter wool
- Contents: 1 pack filter wool for aquarium filters



### You may also be interested in

You can find a complete overview here: https://www.jbl.de/qr/62311



JBL FilterPad VL Cotton fleece for CristalProfi aquarium filters



JBL FilterPad F15 Coarse foam pad for aquarium filter CristalProfi



JBL FilterPad F35 Fine foam pad for aquarium filter CristalProfi




JBL CombiBloc CristalProfi e Pre-filter pads and filter foam for CristalProfi e









Filter materials serve multiple purposes:

They either filter suspended matter out of the water (mechanical filtering), provide settlement space for bacteria that break down pollutants (biological filtering), or remove harmful substances from the water (shamist leber materials)

(chemical/absorptive filter materials). The filter material should always be fitted from coarse to fine inside the filter, as fine filter material clogs too quickly at the beginning. A coarse filter sponge or filter ceramic (JBL Cermec) is ideal for the first filter stage. For biological filtering, JBL Sintomec and JBL Micromec provide the largest possible surface area for pollutant-degrading bacteria. Activated carbon (JBL Carbomec active and Carbomec ultra) filters out discolouration and medication residues. JBL NitratEx and JBL BioNitrat Ex filter out algae promoting nitrates. JBL PhosEx Ultra helps against algae-promoting phosphate. Silicat Ex ultra helps against silicic acid which promotes diatoms.

Filter floss should be used as the last filter material in external filters so that it does not get dirty too quickly. Only in multi-chamber biofilters can filter floss also be used as the first medium, as it is visible there and can be replaced in a few minutes. Filter floss cannot simply be washed out and should always be replaced during filter cleaning.

Absolutely safe

No pollutants released into the aquarium due to residue-free floss.

Further information	
FAQ	~
Blog	~
Press	~
Laboratory/calculator	×
Worth reading	~
Spare parts	×
Video	~
GarantiePlus	×
Instructions	~
QR code	







# JBL Symec Filter floss Product details




Article data					
Product name	JBL Symec Filter Floss 100 g	JBL Symec Filter Floss 250 g	JBL Symec Filter Floss 500 g	JBL Symec Filter Floss 1000 g	
Art. No.	6231100	6231300	6231500	6231700	
EAN number	4014162623119	4014162623133	4014162623157	4014162623171	
EAN as barcode					
Content	100 g	250 g	500 g	1000 g	
-	-	-	-	-	
Expiry months	-	-	-	-	
RRP incl. VAT	3,03 €	6,33 €	12,10 €	20,50 €	
Base price	30.3 €	25.32 €	24.2 €	20.5 €	
Nominal filling quantity	100 g	250 g	500 g	1000 g	
Base quantity	1 kg	1 kg	1 kg	1 kg	
Gross weight	103 g	256 g	510 g	1030 g	
Net weight	100 g	250 g	500 g	1000 g	
Weight change	1000	1000	1000	1000	

Disposal				
Product name	JBL Symec Filter Floss 100 g	JBL Symec Filter Floss 250 g	JBL Symec Filter Floss 500 g	JBL Symec Filter Floss 1000 g
Art. No.	6231100	6231300	6231500	6231700
Green dot	✓	✓	✓	✓
Group electronic waste	-	-	-	-
Disposal weight	-	-	-	-
Battery type	-	-	-	-
Battery return	-	-	-	-
Battery rechargeable	-	-	-	-
Disposal weight battery	-	-	-	-
Non-returnable glass	-	-	-	-
PPK	-	-	-	-
Plastic small	8 g	15 g	25 g	18 g
Plastic large	-	-	-	-
Disposal weight metal	0 g	0 g	0 g	0 g

Features				
Product name	JBL Symec Filter Floss 100 g	JBL Symec Filter Floss 250 g	JBL Symec Filter Floss 500 g	JBL Symec Filter Floss 1000 g
Art. No.	6231100	6231300	6231500	6231700
Animal species	Arowana, Axolotl, Barbels, Bettas, Bichirs/reedfish, Blowfish, Catfish, Cichlids (South America), Corals, Crayfish, Danions, Discus, Dwarf shrimps, Flowerhorn, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, freshwater butterflyfish, turtle	Arowana, Axolotl, Barbels, Bettas, Bichirs/reedfish, Blowfish, Catfish, Cichlids (South America), Corals, Crayfish, Danions, Discus, Dwarf shrimps, Flowerhorn, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, freshwater butterflyfish, turtle	Arowana, Axolotl, Barbels, Bettas, Bichirs/reedfish, Blowfish, Catfish, Cichlids (South America), Corals, Crayfish, Danions, Discus, Dwarf shrimps, Flowerhorn, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, freshwater butterflyfish, turtle	Arowana, Axolotl, Barbels, Bettas, Bichirs/reedfish, Blowfish, Catfish, Cichlids (South America), Corals, Crayfish, Danions, Discus, Dwarf shrimps, Flowerhorn, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, freshwater butterflyfish, turtle

Date: 31.12.2023 Produced by:








## TECHNICAL PRODUCTS TECHNIK TECHNIQUE





Features				
Animal size	For all animal sizes			
Animal age group	All aquarium fish	All aquarium fish	All aquarium fish	All aquarium fish
Volume of habitat	For all aquariums	For all aquariums	For all aquariums	For all aquariums
Material	Polyester	Polyester	Polyester	Polyester
Food type	-	-	-	-
Colour	white	white	white	white
Dosage	-	-	-	-
Transport conditions	-	-	-	-

Date: 31.12.2023 Produced by:








Electronic label / illuminant				
Product name	JBL Symec Filter Floss 100 g	JBL Symec Filter Floss 250 g	JBL Symec Filter Floss 500 g	JBL Symec Filter Floss 1000 g
Art. No.	6231100	6231300	6231500	6231700
Ambient temperature	-	-	-	-
Start time	-	-	-	-
Mercury	-	-	-	-
Tube length	-	-	-	-
Service life	-	-	-	-
Lumen	-	-	-	-
CRI value	-	-	-	-
Dimmable	-	-	-	-
Switching cycles	-	-	-	-
PAR value	-	-	-	-
Energy efficiency class	-	-	-	-
UV-A	-	-	-	-
UV-B	-	-	-	-
UV-C	-	-	-	-
Colour temperature	-	-	-	-
Base designation	-	-	-	-

Technical data				
Product name	JBL Symec Filter Floss 100 g	JBL Symec Filter Floss 250 g	JBL Symec Filter Floss 500 g	JBL Symec Filter Floss 1000 g
Art. No.	6231100	6231300	6231500	6231700
Range in litres	-	-	-	-
Range from - to	-	-	-	-
Range in days	-	-	-	-
Range tank length	-	-	-	-
Output in watts	-	-	-	-
Output per hour	-	-	-	-
Output per day	-	-	-	-
Height	285 mm	350 mm	510 mm	630 mm
Length	180 mm	200 mm	205 mm	280 mm
Width	80 mm	100 mm	105 mm	130 mm
Diameter	-	-	-	-
Voltage	-	-	-	-
For	-	-	-	-
T8 26mm (watt)	-	-	-	-
T5 16mm (watt)	-	-	-	-
Size	-	-	-	-
Content for	-	-	-	-
Filter container volume	-	-	-	-
Volume filter media	-	-	-	-
Hose connections pressure/out	-	-	-	-
Hose connections suction/in	-	-	-	-
Delivery head	-	-	-	-









Food type	-
Sub product type	-
Dosing	-



Date: 31.12.2023 Produced by:







## Additional information for the specialist trade sector

Article data				
Product name	JBL Symec Filter Floss 100 g	JBL Symec Filter Floss 250 g	JBL Symec Filter Floss 500 g	JBL Symec Filter Floss 1000 g
Art. No.	6231100	6231300	6231500	6231700
VAT	19%	19%	19%	19%
Sales unit (SU)	1	1	1	1
Volume packaging	4.11	71	111	231
Dimensions (l/w/h)	80 mm/180 mm/285 mm	100 mm/200 mm/350 mm	105 mm/205 mm/510 mm	130 mm/280 mm/630 mm
Layer	100	0	0	0
Pallet	400	280	168	84
Category of products	2	2	2	2
Customs tariff	56012290	56012290	56012290	56012290
Country of origin	DE	DE	DE	DE
Type of packaging	Bag / sack	Bag / sack	Bag / sack	Bag / sack

PU 1 data				
Product name	JBL Symec Filter Floss 100 g	JBL Symec Filter Floss 250 g	JBL Symec Filter Floss 500 g	JBL Symec Filter Floss 1000 g
Art. No.	6231100	6231300	6231500	6231700
PU 1 material	cardboard gr	cardboard gr	cardboard gr	cardboard gr
PU 1 weight	1070 g	1210 g	1210 g	1210 g
PU 1 lengh	598 mm	528 mm	528 mm	528 mm
PU 1 width/depth	388 mm	388 mm	388 mm	388 mm
PU 1 height	455 mm	635 mm	635 mm	635 mm

PU 2 data				
Product name	JBL Symec Filter Floss 100 g	JBL Symec Filter Floss 250 g	JBL Symec Filter Floss 500 g	JBL Symec Filter Floss 1000 g
Art. No.	6231100	6231300	6231500	6231700
PU 2 material	None	None	None	None
PU 2 weight	-	-	-	-
PU 2 lengh	-	-	-	-
PU 2 width/depth	-	-	-	-
PU 2 height	-	-	-	-

Trade data				
Product name	JBL Symec Filter Floss 100 g	JBL Symec Filter Floss 250 g	JBL Symec Filter Floss 500 g	JBL Symec Filter Floss 1000 g
Art. No.	6231100	6231300	6231500	6231700
Till receipt text	JBL Symec 100g	JBL Symec 250g	JBL Symec 500g	JBL Symec 1000g
Shelf placement	-	-	-	-

31.12.2023 Date: Produced by:

Country: Gibraltar Language: English

