



Quick test to determine the total hardness







- Simple and reliable control of water values in ponds/aquariums. Determines the ideal total hardness for freshwater aquariums and ponds
- Colour change test: fill plastic cuvette with sample water, add reagent drop by drop until colour changes from red to green. Number of drops = level of GH (total hardness)
- When to use: for setting up a new aquarium: once a week
- Online Laboratory: regular control for a healthy aquarium/pond with conditions close to nature. JBL has water tests in the form of quick tests or colour change tests for every water analysis. Testing aquarium water for healthy, clear water
- Package contents: 1 quick test, GH Test. Incl. reagent and plastic vial. Refill reagent avaiable separately



You may also be interested in

You can find a complete overview here: https://www.jbl.de/gr/25350





JBL Test Combi Set plus Fe

Test case with most important water tests + iron test



JBL Testlab Test case with 13 tests f. freshwater analysis Testlab



JBL Testlab Marin Professional test case for the analysis of saltwater



JBL EasyTest 6in1 Test strips for quick aquarium water testing



Date: 31.12.2023

Produced by:











JBL component set for water tests



Country: Georgia

Language: English











JBL assortment Test Colour Charts

Colour charts for JBL water tests

JBL Testlab empty + inserts



Country: Georgia

Language: English







Product information

Healthy aquarium/healthy pond with conditions close to nature The right water values are dependent on the fish stock and the plants in the aquarium/pond. Even if the water looks clear it can be contaminated. With bad values diseases or algae can appear in the aquarium/pond. To maintain a healthy aquarium/pond with conditions close to nature it is important to check and adapt the water values regularly.

For each water analysis JBL provides water tests in the form of quick tests or colour change tests. These determine a certain value or several values in one go. With these water tests you can recognise algae problems and negative nitrate, nitrite, potassium, magnesium values etc.

Depending on origin and consistency of the subsoil, water may contain varying quantities of alkaline earth salts, usually in the form of calcium and magnesium salts. General hardness, by definition, is the sum of calcium and magnesium ions suspended in the water. Most fish and plants thrive well at a general hardness of between roughly 8 and 20°dH. Excessive total hardness values, should, however, be counteracted by the appropriate measures.

JBL Online Laboratory

Regular control for a healthy aquarium/pond with conditions close to nature. Enter your water values into the JBL Online Laboratory and get an in-depth analysis of your values within seconds.

Recommended general hardness values:

Freshwater aquarium (community aquarium): 8-25 °dGH
Lake Malawi/Lake Tanganyika aquarium: 5-20 °dGH
Plant aquarium with few fish (aquascaping): 3-10 °dGH
Marine aquarium: not measurable with GH test because GH too high. Please test calcium and

magnesium separately

Pond: 6-20 °dGH

| Further information | |
|-----------------------|---|
| FAQ | ~ |
| Blog | ~ |
| Press | ~ |
| Laboratory/calculator | × |
| Worth reading | ~ |
| Spare parts | ~ |
| Video | × |
| GarantiePlus | × |
| Instructions | ~ |
| QR code | |









Product details

| Article data | | |
|--------------------------|-----------------|----------------|
| Product name | JBL GH Test Set | JBL GH Reagent |
| Art. No. | 2535000 | 2535100 |
| EAN number | 4014162253507 | 4014162253514 |
| EAN as barcode | | |
| Content | - | - |
| - | - | - |
| Expiry months | 24 | 24 |
| RRP incl. VAT | - | - |
| Base price | - | - |
| Nominal filling quantity | - | - |
| Base quantity | 1 | 1 |
| Gross weight | 36 g | 27 g |
| Net weight | 24 g | 10 g |
| Weight change | 1000 | 1000 |

| Disposal | | |
|-------------------------|-----------------|----------------|
| Product name | JBL GH Test Set | JBL GH Reagent |
| Art. No. | 2535000 | 2535100 |
| Green dot | ✓ | ✓ |
| Group electronic waste | - | - |
| Disposal weight | - | - |
| Battery type | - | - |
| Battery return | - | - |
| Battery rechargeable | - | - |
| Disposal weight battery | - | - |
| Non-returnable glass | - | - |
| PPK | 14 g | 7 g |
| Plastic small | 6 g | - |
| Plastic large | - | - |
| Disposal weight metal | 0 g | 0 g |

| Features | | |
|----------------|---|---|
| Product name | JBL GH Test Set | JBL GH Reagent |
| Art. No. | 2535000 | 2535100 |
| Animal species | Arowana, Axolotl, Barbels, Bettas, Bichirs/reedfish, Blowfish, Catfish, Cichlids (South America), Crayfish, Danions, Discus, Dwarf shrimps, Flowerhorn, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, freshwater butterflyfish, turtle | Arowana, Axolotl, Barbels, Bettas, Bichirs/reedfish, Blowfish, Catfish, Cichlids (South America), Crayfish, Danions, Discus, Dwarf shrimps, Flowerhorn, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, freshwater butterflyfish, turtle |



31.12.2023 Date: Produced by:

Country: Georgia

Language: English





| Features | | |
|----------------------|--|--|
| Animal size | For all animal sizes | For all animal sizes |
| Animal age group | All aquarium fish | All aquarium fish |
| Volume of habitat | For all aquariums | For all aquariums |
| Material | Ammonium hydroxide, ethylenediaminetetraacetic acid (free acid), eriochrome black, purified water, ethanol | Ammonium hydroxide, ethylenediaminetetraacetic acid (free acid), eriochrome black, purified water, ethanol |
| Food type | - | - |
| Colour | red | red |
| Dosage | 5 ml sample water, add reagent drop by drop, count drops, sway after each drop until colour change takes place from red to green | 5 ml sample water, add reagent drop by drop, count drops, sway after each drop until colour change takes place from red to green |
| Transport conditions | - | - |



Country: Georgia Language: English





| Electronic label / illuminant | | |
|-------------------------------|-----------------|----------------|
| Product name | JBL GH Test Set | JBL GH Reagent |
| Art. No. | 2535000 | 2535100 |
| Ambient temperature | - | - |
| Start time | - | - |
| Mercury | - | - |
| Tube length | - | - |
| Service life | - | - |
| Lumen | - | - |
| CRI value | - | - |
| Dimmable | - | - |
| Switching cycles | - | - |
| PAR value | - | - |
| Energy efficiency class | - | - |
| UV-A | - | - |
| UV-B | - | - |
| UV-C | - | - |
| Colour temperature | - | - |
| Base designation | - | - |

| Technical data | | |
|-------------------------------|-----------------|----------------|
| Product name | JBL GH Test Set | JBL GH Reagent |
| Art. No. | 2535000 | 2535100 |
| Range in litres | - | - |
| Range from - to | - | - |
| Range in days | - | - |
| Range tank length | - | - |
| Output in watts | - | - |
| Output per hour | - | - |
| Output per day | - | - |
| Height | - | - |
| Length | - | - |
| Width | - | - |
| Diameter | - | - |
| Voltage | - | - |
| For | - | - |
| T8 26mm (watt) | - | - |
| T5 16mm (watt) | - | - |
| Size | - | - |
| Content for | - | - |
| Filter container volume | - | - |
| Volume filter media | - | - |
| Hose connections pressure/out | - | - |
| Hose connections suction/in | - | - |
| Delivery head | - | - |



Date: 31.12.2023 Produced by:







| Food type | - |
|------------------|--|
| Sub product type | - |
| Dosing | 5 ml sample water, add reagent drop by drop, count drops, sway after each drop until colour change takes place from red to green |











A Safety information

| - safety information in accordance with GHS | | |
|---|--|--|
| Signal word | Danger | |
| Hazard symbol | | |
| Risk phrase | H226: Flammable liquid and vapour. H314: Causes severe skin burns and eye damage. | |
| safety note | P102: Keep out of reach of children. P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P261: Avoid breathing dust/fume/gas/mist/vapours/spray. P280: Wear protective gloves/protective clothing/eye protection/face protection. P303+P361+P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P310: Immediately call a POISON CENTER/doctor. | |

| Further safety information | |
|----------------------------|---|
| Safety instructions | ✓ |
| QR code | |



31.12.2023 Date: Produced by:







Additional information for the specialist trade sector

| Article data | | |
|----------------------|----------------------|-------------------------------|
| Product name | JBL GH Test Set | JBL GH Reagent |
| Art. No. | 2535000 | 2535100 |
| VAT | 19% | 19% |
| Sales unit (SU) | 6 | 1 |
| Volume packaging | 0.3I | 0.451 |
| Dimensions (l/w/h) | 35 mm/70 mm/125 mm | 30 mm/107 mm/140 mm |
| Layer | 960 | 0 |
| Pallet | 2880 | 0 |
| Category of products | 1 | 1 |
| Customs tariff | 38220000 | 38220000 |
| Country of origin | DE | DE |
| Type of packaging | Faltschachtel/Karton | Skin / blister card packaging |

| PU 1 data | | |
|------------------|-----------------|----------------|
| Product name | JBL GH Test Set | JBL GH Reagent |
| Art. No. | 2535000 | 2535100 |
| PU 1 material | film gr | cardboard gr |
| PU 1 weight | 4.3 g | 660 g |
| PU 1 lengh | 215 mm | 580 mm |
| PU 1 width/depth | 35 mm | 390 mm |
| PU 1 height | 185 mm | 310 mm |

| PU 2 data | | |
|------------------|-----------------|----------------|
| Product name | JBL GH Test Set | JBL GH Reagent |
| Art. No. | 2535000 | 2535100 |
| PU 2 material | cardboard gr | None |
| PU 2 weight | 500 g | - |
| PU 2 lengh | 572 mm | - |
| PU 2 width/depth | 377 mm | - |
| PU 2 height | 410 mm | - |

| Trade data | | |
|-------------------|-----------------|----------------|
| Product name | JBL GH Test Set | JBL GH Reagent |
| Art. No. | 2535000 | 2535100 |
| Till receipt text | GH Test | GH Test Re |
| Shelf placement | - | - |



31.12.2023 Date: Produced by: