



JBL GoldPearls mini CLICK

Premium quality for small veiltails

Suitable for: 



- Complete nutrition, ideal growth and vibrant colours for small goldfish and veiltails
- Nutritious and easily to digest: fine food pellets create ideal feeding conditions for small, low swimming fish (min. 4 cm). No air swallowing thanks to sinking food pellets.
- No water clouding: reduced algae growth thanks to balanced phosphate content, better water quality due to better digestibility of the food, which reduces fish excrement
- Fish choose JBL food: 98.5% of all fish species ate the JBL food immediately during research expeditions in freshwaters. No processing of cheap fish meal, fish meat from fish fillet production for humans is used
- Contents: premium main food, with the Click doser for aquarium fish. After opening to be used up within 4 months = 1 click feeds 5 fish. In high-quality packaging with foil seal, unopened can be stored for 3 years

 You may also be interested in

You can find a complete overview here: <https://www.jbl.de/qr/40643>



JBL Holiday Red
Complete holiday food for goldfish





JBL GoldPearls mini CLICK

Accessories



JBL Atvitol
Multivitamin drops for
aquarium fish



JBL AutoFood BLACK
Black automatic feeder for
aquarium fish



JBL AutoFood WHITE
White automatic feeder for
aquarium fish



JBL GoldPearls mini CLICK



Product information

Animal welfare is important to us

No processing of cheap fish meal, uses fish meat from fish fillet production with the motto: The large fillet for humans and the small fillet for our aquarium fish.

For years JBL has supported the largest international organisation for the protection of sharks "Shark Project".

Other good reasons to use JBL fish food:

- The use of pure fish protein without cheap fish meal
- Optimal protein/fat ratio
- Mainly proteins from water animals
- Reduction of algae growth and optimum fish growth thanks to adapted phosphate content
- Impressive acceptance on the part of fish: research expeditions with fish feeding trials in the wild
- Very slight vitamin loss due to airtight sealed packaging


Ahead through research

The results of the JBL research expeditions, combined with the expertise of the JBL research and development team, has resulted in optimal and balanced food mixtures made of high-quality ingredients.

Feeding recommendation

Feed 1 – 2 times a day, as much as can be eaten within a few minutes. Young, growing fish 3 – 4 times a day in the same way.

Unopened the food tins can be kept for 3 years, after opening use up within 4 months, since valuable vitamins will deteriorate. Choose the size of the tin (100 ml, 250 ml) in accordance with the monthly consumption stated above and the number of the fish in the aquarium. Complements complete food and available separately: complementary food and treats.

| Further information | |
|-----------------------|---|
| FAQ | ✓ |
| Blog | ✓ |
| Press | ✓ |
| Laboratory/calculator | ✗ |
| Worth reading | ✓ |
| Spare parts | ✗ |
| Video | ✓ |
| GarantiePlus | ✗ |
| Instructions | ✓ |
| QR code |  |






JBL GoldPearls mini CLICK



Product details

| Article data | |
|--------------------------|---|
| Product name | JBL GoldPearls mini 100 ml CLICK |
| Art. No. | 4064300 |
| EAN number | 4014162406439 |
| EAN as barcode |  |
| Weight | 56 g |
| Size | 100 ml |
| Expiry months | 36 |
| RRP incl. VAT | 6,20 € |
| Base price | 110.71 € |
| Nominal filling quantity | 56 g |
| Base quantity | 1 kg |
| Gross weight | 72.6 g |
| Net weight | 56 g |
| Weight change | 1000 |

| Disposal | |
|-------------------------|----------------------------------|
| Product name | JBL GoldPearls mini 100 ml CLICK |
| Art. No. | 4064300 |
| Green dot | ✓ |
| Group electronic waste | - |
| Disposal weight | - |
| Battery type | - |
| Battery return | - |
| Battery rechargeable | - |
| Disposal weight battery | - |
| Non-returnable glass | - |
| PPK | - |
| Plastic small | 25 g |
| Plastic large | - |
| Disposal weight metal | 0 g |

| Features | |
|-------------------|---|
| Product name | JBL GoldPearls mini 100 ml CLICK |
| Art. No. | 4064300 |
| Animal species | Veiltails |
| Animal size | min.4 cm |
| Animal age group | Small veiltails |
| Volume of habitat | For all aquariums |
| Material | You will find details of the materials used under Details & Composition |





| Features | |
|----------------------|---|
| Food type | granulates |
| Colour | brown |
| Dosage | Feed 1 – 2 times a day, as much as can be eaten within a few minutes. Young, growing fish 3 – 4 times a day in the same way |
| Transport conditions | - |





JBL GoldPearls mini CLICK

| Electronic label / illuminant | |
|-------------------------------|----------------------------------|
| Product name | JBL GoldPearls mini 100 ml CLICK |
| Art. No. | 4064300 |
| Ambient temperature | - |
| Start time | - |
| Mercury | - |
| Tube length | - |
| Service life | - |
| Lumen | - |
| CRI value | - |
| Dimmable | - |
| Switching cycles | - |
| PAR value | - |
| Energy efficiency class | - |
| UV-A | - |
| UV-B | - |
| UV-C | - |
| Colour temperature | - |
| Base designation | - |

| Technical data | |
|-------------------------------|----------------------------------|
| Product name | JBL GoldPearls mini 100 ml CLICK |
| Art. No. | 4064300 |
| Range in litres | - |
| Range from - to | - |
| Range in days | - |
| Range tank length | - |
| Output in watts | - |
| Output per hour | - |
| Output per day | - |
| Height | - |
| Length | - |
| Width | - |
| Diameter | - |
| Voltage | - |
| For | - |
| T8 26mm (watt) | - |
| T5 16mm (watt) | - |
| Size | - |
| Content for | - |
| Filter container volume | - |
| Volume filter media | - |
| Hose connections pressure/out | - |
| Hose connections suction/in | - |
| Delivery head | - |






JBL GoldPearls mini CLICK


| | |
|------------------|---|
| Food type | granulates |
| Sub product type | complete food for ornamental fish |
| Dosing | Feed 1 - 2 times a day, as much as can be eaten within a few minutes. Young, growing fish 3 - 4 times a day in the same way |

Analytical components



| | |
|---------------|--------|
| Crude protein | 38 % |
| Fat content | 4 % |
| Crude fibre | 4.5 % |
| Crude ash | 10.5 % |

Composition



| |
|---------------------------|
| Fish and fish by-products |
| Molluscs and crustaceans |
| Vegetable by-products |
| Vegetables |
| Cereals |
| Algae |
| Yeast |

Additives

| |
|----------------------------------|
| Colourings |
| Astaxanthin E161j food colouring |

| |
|------------------------------------|
| Antioxidants |
| E 306 (natural vitamin E extracts) |

| | |
|---|-------------|
| Vitamins, provitamins and other chemically defined substances having a similar effect (per 1000 g) | |
| Vitamin A | 25000 I. E. |
| Vitamin D ₃ | 2000 I. E. |
| Vitamin E | 300 mg |
| Vitamin C(stable) | 200 mg |





Additional information for the specialist trade sector

| Article data | |
|----------------------|----------------------------------|
| Product name | JBL GoldPearls mini 100 ml CLICK |
| Art. No. | 4064300 |
| VAT | 7% |
| Sales unit (SU) | 6 |
| Volume packaging | 0.26l |
| Dimensions (l/w/h) | 50 mm/50 mm/105 mm |
| Layer | 672 |
| Pallet | 2688 |
| Category of products | 1 |
| Customs tariff | 23099041 |
| Country of origin | DE |
| Type of packaging | can |

| PU 1 data | |
|------------------|----------------------------------|
| Product name | JBL GoldPearls mini 100 ml CLICK |
| Art. No. | 4064300 |
| PU 1 material | film gr |
| PU 1 weight | 3.6 g |
| PU 1 length | 92 mm |
| PU 1 width/depth | 142 mm |
| PU 1 height | 105 mm |

| PU 2 data | |
|------------------|----------------------------------|
| Product name | JBL GoldPearls mini 100 ml CLICK |
| Art. No. | 4064300 |
| PU 2 material | cardboard gr |
| PU 2 weight | 320 g |
| PU 2 length | 390 mm |
| PU 2 width/depth | 295 mm |
| PU 2 height | 175 mm |

| Trade data | |
|-------------------|----------------------------------|
| Product name | JBL GoldPearls mini 100 ml CLICK |
| Art. No. | 4064300 |
| Till receipt text | GoldPearls mini |
| Shelf placement | - |