




# JBL PROSCAN RECHARGE

Refill test strips for analysis via smartphone

Suitable for:



- Water analysis and diagnosis via smartphone: easy and safe control of the water values in freshwater aquariums. For the measurement of: general hardness, carbonate hardness, pH value, nitrite, nitrate, chlorine and CO<sub>2</sub> calculation
- The new generation of water analysis: download ProScan app free of charge, insert test strip into aquarium water, position strip on colour chart for analysis, scan colour chart, values are determined
- Fast diagnosis – precise results: test value results are displayed as numbers and additionally evaluated (good/medium/poor). Additional recommendations for optimal values in freshwater aquariums and ponds
- Compatible with iPhone and iPad: iOS 5.0 or higher, Android 4.0 or higher and a camera with autofocus needed.
- Contents: 1 pack ProScan Recharge with 24 water analysis strips



You may also be interested in

You can find a complete overview here: <https://www.jbl.de/qr/25421>



## JBL PROAQUATEST COMBISET Plus Fe

Test case with most important  
water tests + iron test



## JBL PROAQUATEST COMBISET Plus NH<sub>4</sub>

Test case with most important  
water tests + NH<sub>4</sub> test



## JBL PROAQUATEST LAB

Test case with 13 water tests  
for freshwater analysis



## JBL PROAQUATEST LAB PROSCAPE

Test case for water analysis in  
plant aquariums



# JBL PROSCAN RECHARGE



Spare parts



JBL ProScan ColorCard



# JBL PROSCAN RECHARGE



## Product information

Healthy aquariums with conditions modelled on nature

The right water values are dependent on the fish stock and plants in the aquarium. Even if the water looks clear it can be contaminated. With bad values diseases or algae can appear in the aquarium. To maintain a healthy aquarium with conditions close to nature it is important to check and adapt the water values regularly.

Water analysis using your smartphone – the new generation

Fast diagnosis, precise results and recommended solutions

Saves measurement values

Displays measurement evaluation


New developments - new functions and uses

Fast diagnosis - precise results

Download the JBL app free of charge, follow the app's instructions, immerse the test strip into the aquarium water, position the test strip on the JBL colour chart for the analysis, scan colour chart, values will be determined, individual solution recommendations and tips will be transmitted.

Note: The test strips are only suitable for one measurement each. A repeated scan of the strip leads to different values, since the strip goes on reacting indefinitely.

During the analysis "aquarium freshwater" and "water" the current CO<sub>2</sub> content in the water will be calculated as an additional bonus. Depending on the analysis results in the sections "Aquarium Freshwater" and "Garden Pond" individualized recommendations and hints will be transmitted as to which JBL products will help to set the ideal/optimal water values.


Further information	
FAQ	✓
Blog	✓
Press	✓
Laboratory/calculator	✗
Worth reading	✓
Spare parts	✓
Video	✓
GarantiePlus	✗
Instructions	✓
QR code	



# JBL PROSCAN RECHARGE



## Product details

Article data	
Product name	JBL PROSCAN RECHARGE
Art. No.	2542100
EAN number	4014162254214
EAN as barcode	
Content	24 Tests
-	-
Expiry months	30
RRP incl. VAT	22,87 €
Base price	-
Nominal filling quantity	-
Base quantity	1
Gross weight	57 g
Net weight	9 g
Weight change	1000

Disposal	
Product name	JBL PROSCAN RECHARGE
Art. No.	2542100
Green dot	✓
Group electronic waste	-
Disposal weight	-
Battery type	-
Battery return	-
Battery rechargeable	-
Disposal weight battery	-
Non-returnable glass	-
PPK	15 g
Plastic small	15 g
Plastic large	-
Disposal weight metal	0 g

Features	
Product name	JBL PROSCAN RECHARGE
Art. No.	2542100



Features	
Animal species	Armored catfish, Arowana, Axolotl, Barbels, Bettas, Bichirs/ reedfish, Blowfish, Catfish, Cichlids (South America), Clawed frogs, Clawed shrimps, Crayfish, Crustaceans, Danions, Discus, Dwarf crayfish, Dwarf shrimps, Flowerhorn, Gill maggots, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, blood parrot cichlids, freshwater butterflyfish, turtle
Animal size	For all animal sizes
Animal age group	All aquarium fish
Volume of habitat	24 tests
Material	-
Food type	-
Colour	multicoloured
Dosage	Press the "start" button and immerse the analysis strip into the water for 2 to 3 seconds. Shake sideways to remove any adhering water. Position the analysis strip on the ColorChart and wait until the time displayed has expired. Using your smartphone compare the test strip with the complete ColorChart to receive the analysis results.
Transport conditions	-



# JBL PROSCAN RECHARGE

Electronic label / illuminant	
Product name	JBL PROSCAN RECHARGE
Art. No.	2542100
Ambient temperature	-
Start time	-
Mercury	-
Tube length	-
Service life	-
Lumen	-
CRI value	-
Dimmable	-
Switching cycles	-
PAR value	-
Energy efficiency class	-
UV-A	-
UV-B	-
UV-C	-
Colour temperature	-
Base designation	-

Technical data	
Product name	JBL PROSCAN RECHARGE
Art. No.	2542100
Range in litres	-
Range from - to	-
Range in days	-
Range tank length	-
Output in watts	-
Output per hour	-
Output per day	-
Height	-
Length	-
Width	-
Diameter	-
Voltage	-
For	-
T8 26mm (watt)	-
T5 16mm (watt)	-
Size	-
Content for	-
Filter container volume	-
Volume filter media	-
Hose connections pressure/out	-
Hose connections suction/in	-
Delivery head	-



# JBL PROSCAN RECHARGE

Food type	-
Sub product type	-
Dosing	Press the "start" button and immerse the analysis strip into the water for 2 to 3 seconds. Shake sideways to remove any adhering water. Position the analysis strip on the ColorChart and wait until the time displayed has expired. Using your smartphone compare the test strip with the complete ColorChart to receive the analysis results.



## Additional information for the specialist trade sector

Article data	
Product name	JBL PROSCAN RECHARGE
Art. No.	2542100
VAT	19%
Sales unit (SU)	6
Volume packaging	0.37l
Dimensions (l/w/h)	50 mm/50 mm/150 mm
Layer	540
Pallet	2700
Category of products	1
Customs tariff	38221100
Country of origin	DE
Type of packaging	Faltschachtel/Karton

PU 1 data	
Product name	JBL PROSCAN RECHARGE
Art. No.	2542100
PU 1 material	film gr
PU 1 weight	40 g
PU 1 length	50 mm
PU 1 width/depth	300 mm
PU 1 height	150 mm

PU 2 data	
Product name	JBL PROSCAN RECHARGE
Art. No.	2542100
PU 2 material	cardboard gr
PU 2 weight	520 g
PU 2 length	460 mm
PU 2 width/depth	315 mm
PU 2 height	325 mm

Trade data	
Product name	JBL PROSCAN RECHARGE
Art. No.	2542100
Till receipt text	ProScan Recharge
Shelf placement	-