




JBL PROCRISTAL i30 FilterSponge

Replacement filter sponge for CristalProfi i 30

Suitable for:



- Clear water, healthy fish: replacement filter sponge for JBL ProCristal i 30
- Clear aquarium water: replace the filter sponge every 3 months.
- Best cleaning: pore size 30 ppi
- Best biological filtering: promotes the settlement of beneficial cleansing bacteria.
- Package contents: replacement sponge for JBL CristalProfi i 30, ProCristal i30 FilterSponge



You may also be interested in

You can find a complete overview here: <https://www.jbl.de/qr/60993>



JBL PROCRISTAL i30 FilterModul

Filter module for ProCristal i30



JBL PROCRISTAL i30 FilterSponge



Product information

Clear aquarium water

Replacing the filter sponge regularly ensures clear aquarium water. The water is perfectly cleaned thanks to the pore size of 30 ppi. At the same time the sponge promotes the settlement of beneficial cleansing bacteria for an ideal biological filtering.

Quality of foam:

Foams used in aquaristics with live animals must always be subjected to a pollutant test, as most plastics release problematic plasticisers into the water. All the plastics used by JBL have passed this test and do NOT release any harmful substances into the water.


Cleaning of plastics:

For many months it is sufficient to rinse the dirty filter sponges under lukewarm running water and then inoculate them with filter bacteria (JBL FilterStart) before using them again. However, as mechanical cleaning NEVER leads to complete cleaning and dirt will remain in the filter sponge, it is recommended you replace the filter sponges about every 6 months.

Inoculating filter material with cleansing bacteria:

Initial colonisation of filter material with cleansing bacteria takes about 1-2 weeks in freshwater and twice as long in seawater. If the filter material was already in use, the time for a "recolonisation" is doubled, as recolonisation does not take place so easily on dead bacterial lawns. The addition of concentrated cleaning bacteria such as in JBL FilterStart shortens the biological start-up time to about 24 hours.

Another reason is to fill the existing settlement space in the filter material with beneficial bacteria so that pathogenic germs find no place to spread. This has now been scientifically proven.

| Further information | |
|-----------------------|---|
| FAQ | ✓ |
| Blog | ✓ |
| Press | ✓ |
| Laboratory/calculator | ✗ |
| Worth reading | ✓ |
| Spare parts | ✗ |
| Video | ✓ |
| GarantiePlus | ✗ |
| Instructions | ✓ |
| QR code |  |



JBL PROCRISTAL i30 FilterSponge



Product details

| Article data | |
|--------------------------|---------------------------------|
| Product name | JBL PROCRISTAL i30 FilterSponge |
| Art. No. | 6099300 |
| EAN number | 4014162609939 |
| EAN as barcode | |
| For | PC i30 |
| - | - |
| Expiry months | - |
| RRP incl. VAT | 4,99 € |
| Base price | - |
| Nominal filling quantity | - |
| Base quantity | 1 |
| Gross weight | 16 g |
| Net weight | 1 g |
| Weight change | 1000 |

| Disposal | |
|-------------------------|---------------------------------|
| Product name | JBL PROCRISTAL i30 FilterSponge |
| Art. No. | 6099300 |
| Green dot | ✓ |
| Group electronic waste | - |
| Disposal weight | - |
| Battery type | - |
| Battery return | - |
| Battery rechargeable | - |
| Disposal weight battery | - |
| Non-returnable glass | - |
| PPK | 6 g |
| Plastic small | 8.5 g |
| Plastic large | - |
| Disposal weight metal | 0 g |

| Features | |
|----------------|--|
| Product name | JBL PROCRISTAL i30 FilterSponge |
| Art. No. | 6099300 |
| Animal species | Arowana, Axolotl, Barbels, Bettas, Bichirs/reedfish, Blowfish, Catfish, Cichlids (South America), Crayfish, Danions, Discus, Dwarf shrimps, Flowerhorn, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, freshwater butterflyfish, turtle |



| Features | |
|----------------------|------------------------|
| Animal size | 1 - 5 cm |
| Animal age group | Small aquarium fish |
| Volume of habitat | 10 - 40 L (15 - 40 cm) |
| Material | polyether |
| Food type | - |
| Colour | blue |
| Dosage | - |
| Transport conditions | - |



JBL PROCRISTAL i30 FilterSponge

| Electronic label / illuminant | |
|-------------------------------|---------------------------------|
| Product name | JBL PROCRISTAL i30 FilterSponge |
| Art. No. | 6099300 |
| Ambient temperature | - |
| Start time | - |
| Mercury | - |
| Tube length | - |
| Service life | - |
| Lumen | - |
| CRI value | - |
| Dimmable | - |
| Switching cycles | - |
| PAR value | - |
| Energy efficiency class | - |
| UV-A | - |
| UV-B | - |
| UV-C | - |
| Colour temperature | - |
| Base designation | - |

| Technical data | |
|-------------------------------|---------------------------------|
| Product name | JBL PROCRISTAL i30 FilterSponge |
| Art. No. | 6099300 |
| Range in litres | - |
| Range from - to | - |
| Range in days | - |
| Range tank length | - |
| Output in watts | - |
| Output per hour | - |
| Output per day | - |
| Height | - |
| Length | - |
| Width | - |
| Diameter | - |
| Voltage | - |
| For | - |
| T8 26mm (watt) | - |
| T5 16mm (watt) | - |
| Size | - |
| Content for | - |
| Filter container volume | - |
| Volume filter media | - |
| Hose connections pressure/out | - |
| Hose connections suction/in | - |
| Delivery head | - |



JBL PROCRISTAL i30 FilterSponge

| | |
|------------------|---|
| Food type | - |
| Sub product type | - |
| Dosing | - |



Additional information for the specialist trade sector

| Article data | |
|----------------------|---------------------------------|
| Product name | JBL PROCRISTAL i30 FilterSponge |
| Art. No. | 6099300 |
| VAT | 19% |
| Sales unit (SU) | 4 |
| Volume packaging | 0.63l |
| Dimensions (l/w/h) | 50 mm/96 mm/140 mm |
| Layer | 0 |
| Pallet | 0 |
| Category of products | 3 |
| Customs tariff | 84219990 |
| Country of origin | IT |
| Type of packaging | Skin / blister card packaging |

| PU 1 data | |
|------------------|---------------------------------|
| Product name | JBL PROCRISTAL i30 FilterSponge |
| Art. No. | 6099300 |
| PU 1 material | cardboard gr |
| PU 1 weight | 44 g |
| PU 1 length | 120 mm |
| PU 1 width/depth | 100 mm |
| PU 1 height | 145 mm |

| PU 2 data | |
|------------------|---------------------------------|
| Product name | JBL PROCRISTAL i30 FilterSponge |
| Art. No. | 6099300 |
| PU 2 material | cardboard gr |
| PU 2 weight | 642 g |
| PU 2 length | 521 mm |
| PU 2 width/depth | 379 mm |
| PU 2 height | 316 mm |

| Trade data | |
|-------------------|---------------------------------|
| Product name | JBL PROCRISTAL i30 FilterSponge |
| Art. No. | 6099300 |
| Till receipt text | CP i30Sponge |
| Shelf placement | - |