



# JBL PRONOVO BOTIA TAB M

NEW

Food tablets for all loaches from 1-20 cm

Suitable for: 



- Ideal staple food for all loaches (Botia species, sucking loaches, hillstram loaches, Pangio species)
- Use: either let it sink to the bottom of the aquarium or attach it to the aquarium glass with light pressure from the inside
- Natural raw ingredients without artificial additives for a natural diet and healthy aquarium fish. Tested on JBL research expeditions in the tropics
- Very high acceptance even with choosy species thanks to high-quality raw ingredients such as salmon, spirulina and shrimps
- Package contents: original sealed, air- and light-tight, recyclable tin with dosing aid inside the lid (for 100 ml, 250 ml, 1000 ml)



# JBL PRONOVO BOTIA TAB M




## Product information

All loaches come from Asia and are ground-oriented fish. They therefore eat mainly from the substrate. Many species, such as the clown loach, skunk loach, yoyo loach and orange-finned loach are good snail eaters!

How do bottom dwellers in a community aquarium get their food? You'll have seen this: Someone throws a food tablet into the aquarium for their bottom-dwelling loaches, but the fish in the other water zones grab the food tablet and the loaches go empty-handed! The other variant: The loaches are supposed to get some of the other fish's food. But too often so much food is fed that it's enough for the loaches and the algae as well! The loaches may not be able to reach food that lands in the decoration. This can spoil and quickly lead to algae problems. This is how it works better: first feed the fish in the upper water layers and then, when they are in a "feeding frenzy", throw in one or two food tablets for the loaches.

How can you proceed if you keep different bottom dwellers in one aquarium? You can feed 2 types of tablets (e.g. JBL PRONOVO CORYDORAS for armored catfish and JBL PRONOVO BOTIA for loaches) at the same time without any problems. It is even beneficial for the fish, as they are offered more variety and yet each group of fish gets its species-appropriate nutrition.

What is the difference between tablets and wafers? Food tablets (e.g. JBL PRONOVO TAB, JBL PRONOVO BOTIA) are mechanically compressed from fine flake food. Wafers, on the other hand, are thinly cut from a pulp and a "continuous sausage" pressed out of the extruder and dried.

| Further information   |   |
|-----------------------|---|
| FAQ                   | ✓   |
| Blog                  | ✓   |
| Press                 | ✓   |
| Laboratory/calculator | ✗   |
| Worth reading         | ✓   |
| Spare parts           | ✗   |
| Video                 | ✓   |
| GarantiePlus          | ✗   |
| Instructions          | ✓   |
| QR code               |  |





# JBL PRONOVO BOTIA TAB M



## Product details

| Article data             |                             |                             |                              |                              |
|--------------------------|-----------------------------|-----------------------------|------------------------------|------------------------------|
| Product name             | JBL PRONOVO BOTIA TAB 100ml | JBL PRONOVO BOTIA TAB 250ml | JBL PRONOVO BOTIA TAB 1000ml | JBL PRONOVO BOTIA TAB M 5.5l |
| Art. No.                 | 3115273                     | 3115300                     | 3115500                      | 3115600                      |
| EAN number               | 4014162056344               | 4014162311535               | 4014162311559                | 4014162311566                |
| EAN as barcode           |                             |                             |                              |                              |
| Content                  | 100 ml                      | 250 ml                      | 1000 ml                      | 5,5 l                        |
| Weight                   | 58 g                        | 150 g                       | 580 g                        | 2900 g                       |
| Expiry months            | 36                          | 36                          | 36                           | 36                           |
| RRP incl. VAT            | 8,43 €                      | 17,12 €                     | 47,49 €                      | 86,30 €                      |
| Base price               | 145,34 €                    | 114,13 €                    | 81,88 €                      | 29,76 €                      |
| Nominal filling quantity | 58 g                        | 150 g                       | 580 g                        | 2900 g                       |
| Base quantity            | 1 kg                        | 1 kg                        | 1 kg                         | 1 kg                         |
| Gross weight             | 78,94 g                     | 185,9 g                     | 659,956 g                    | 3079 g                       |
| Net weight               | 58 g                        | 150 g                       | 580 g                        | 2900 g                       |
| Weight change            | -                           | -                           | -                            | -                            |

| Disposal                |                             |                             |                              |                              |
|-------------------------|-----------------------------|-----------------------------|------------------------------|------------------------------|
| Product name            | JBL PRONOVO BOTIA TAB 100ml | JBL PRONOVO BOTIA TAB 250ml | JBL PRONOVO BOTIA TAB 1000ml | JBL PRONOVO BOTIA TAB M 5.5l |
| Art. No.                | 3115273                     | 3115300                     | 3115500                      | 3115600                      |
| Green dot               | ✓                           | ✓                           | -                            | ✓                            |
| Group electronic waste  | -                           | -                           | -                            | -                            |
| Disposal weight         | -                           | -                           | -                            | -                            |
| Battery type            | -                           | -                           | -                            | -                            |
| Battery return          | -                           | -                           | -                            | -                            |
| Battery rechargeable    | -                           | -                           | -                            | -                            |
| Disposal weight battery | -                           | -                           | -                            | -                            |
| Non-returnable glass    | -                           | -                           | -                            | -                            |
| PPK                     | -                           | -                           | -                            | -                            |
| Plastic small           | 20,94 g                     | 35,9 g                      | 79,96 g                      | -                            |
| Plastic large           | -                           | -                           | -                            | 179 g                        |
| Disposal weight metal   | 0 g                         | 0 g                         | 0 g                          | 0 g                          |

| Features          |   |   |   |   |
|-------------------|---|---|---|---|
| Product name      | JBL PRONOVO BOTIA TAB 100ml   | JBL PRONOVO BOTIA TAB 250ml   | JBL PRONOVO BOTIA TAB 1000ml  | JBL PRONOVO BOTIA TAB M 5.5l  |
| Art. No.          | 3115273   | 3115300   | 3115500   | 3115600   |
| Animal species    | Loaches   | Loaches   | Loaches   | Loaches   |
| Animal size       | 1 - 20 cm   | 1 - 20 cm   | 1 - 20 cm   | 1 - 20 cm   |
| Animal age group  | All aquarium fish   | All aquarium fish   | All aquarium fish   | All aquarium fish   |
| Volume of habitat | For all aquariums   | For all aquariums   | For all aquariums   | For all aquariums   |
| Material          | You will find details of the materials used under Details & Composition | You will find details of the materials used under Details & Composition | You will find details of the materials used under Details & Composition | You will find details of the materials used under Details & Composition |





| Features             |  |  |  |  |
|----------------------|--|--|--|--|
| Food type            | food tablets   | food tablets   | food tablets   | food tablets   |
| Colour               | multicoloured  | multicoloured  | multicoloured  | multicoloured  |
| Dosage               | Feed 1 – 2 times a day, as much as can be eaten within a few minutes | Feed 1 – 2 times a day, as much as can be eaten within a few minutes | Feed 1 – 2 times a day, as much as can be eaten within a few minutes | Feed 1 – 2 times a day, as much as can be eaten within a few minutes |
| Transport conditions | -  | -  | -  | -  |





# JBL PRONOVO BOTIA TAB M

| Electronic label / illuminant |                             |                             |                              |                              |
|-------------------------------|-----------------------------|-----------------------------|------------------------------|------------------------------|
| Product name                  | JBL PRONOVO BOTIA TAB 100ml | JBL PRONOVO BOTIA TAB 250ml | JBL PRONOVO BOTIA TAB 1000ml | JBL PRONOVO BOTIA TAB M 5.5l |
| Art. No.                      | 3115273                     | 3115300                     | 3115500                      | 3115600                      |
| Ambient temperature           | -                           | -                           | -                            | -                            |
| Start time                    | -                           | -                           | -                            | -                            |
| Mercury                       | -                           | -                           | -                            | -                            |
| Tube length                   | -                           | -                           | -                            | -                            |
| Service life                  | -                           | -                           | -                            | -                            |
| Lumen                         | -                           | -                           | -                            | -                            |
| CRI value                     | -                           | -                           | -                            | -                            |
| Dimmable                      | -                           | -                           | -                            | -                            |
| Switching cycles              | -                           | -                           | -                            | -                            |
| PAR value                     | -                           | -                           | -                            | -                            |
| Energy efficiency class       | -                           | -                           | -                            | -                            |
| UV-A                          | -                           | -                           | -                            | -                            |
| UV-B                          | -                           | -                           | -                            | -                            |
| UV-C                          | -                           | -                           | -                            | -                            |
| Colour temperature            | -                           | -                           | -                            | -                            |
| Base designation              | -                           | -                           | -                            | -                            |

| Technical data                |                             |                             |                              |                              |
|-------------------------------|-----------------------------|-----------------------------|------------------------------|------------------------------|
| Product name                  | JBL PRONOVO BOTIA TAB 100ml | JBL PRONOVO BOTIA TAB 250ml | JBL PRONOVO BOTIA TAB 1000ml | JBL PRONOVO BOTIA TAB M 5.5l |
| Art. No.                      | 3115273                     | 3115300                     | 3115500                      | 3115600                      |
| Range in litres               | -                           | -                           | -                            | -                            |
| Range from - to               | -                           | -                           | -                            | -                            |
| Range in days                 | -                           | -                           | -                            | -                            |
| Range tank length             | -                           | -                           | -                            | -                            |
| Output in watts               | -                           | -                           | -                            | -                            |
| Output per hour               | -                           | -                           | -                            | -                            |
| Output per day                | -                           | -                           | -                            | -                            |
| Height                        | 106 mm                      | 146 mm                      | 198 mm                       | 191 mm                       |
| Length                        | 44 mm                       | 56 mm                       | 90 mm                        | 223 mm                       |
| Width                         | 44 mm                       | 56 mm                       | 90 mm                        | 223 mm                       |
| Diameter                      | 44 mm                       | 56 mm                       | 90 mm                        | 223 mm                       |
| Voltage                       | -                           | -                           | -                            | -                            |
| For                           | -                           | -                           | -                            | -                            |
| T8 26mm (watt)                | -                           | -                           | -                            | -                            |
| T5 16mm (watt)                | -                           | -                           | -                            | -                            |
| Size                          | -                           | -                           | -                            | -                            |
| Content for                   | -                           | -                           | -                            | -                            |
| Filter container volume       | -                           | -                           | -                            | -                            |
| Volume filter media           | -                           | -                           | -                            | -                            |
| Hose connections pressure/out | -                           | -                           | -                            | -                            |
| Hose connections suction/in   | -                           | -                           | -                            | -                            |
| Delivery head                 | -                           | -                           | -                            | -                            |






# JBL PRONOVO BOTIA TAB M

|                  |  |
|------------------|--|
| Food type        | food tablets   |
| Sub product type | complete food for ornamental fish                                    |
| Dosing           | Feed 1 - 2 times a day, as much as can be eaten within a few minutes |

## Analytical components



|               |      |
|---------------|------|
| Crude protein | 33 % |
| Fat content   | 6 %  |
| Crude fibre   | 6 %  |
| Crude ash     | 8 %  |

## Composition



|               |
|---------------|
| spirulina     |
| salmon meal   |
| wheat germs   |
| shrimp meal   |
| rice meal     |
| soy bean meal |
| wheat meal    |
| green meal    |
| Apple fibre   |
| wheat gluten  |
| krill meal    |
| corn meal     |
| yeast extract |
| Potato flakes |
| spinach       |

## Additives

|                      |
|----------------------|
| <b>Colourings</b>    |
| are not added by JBL |

|                     |
|---------------------|
| <b>Antioxidants</b> |
| 1b306(i)            |





## Additional information for the specialist trade sector

| Article data         |                             |                             |                              |                              |
|----------------------|-----------------------------|-----------------------------|------------------------------|------------------------------|
| Product name         | JBL PRONOVO BOTIA TAB 100ml | JBL PRONOVO BOTIA TAB 250ml | JBL PRONOVO BOTIA TAB 1000ml | JBL PRONOVO BOTIA TAB M 5.5l |
| Art. No.             | 3115273                     | 3115300                     | 3115500                      | 3115600                      |
| VAT                  | 7%                          | 7%                          | 7%                           | 7%                           |
| Sales unit (SU)      | 6                           | 6                           | 3                            | 1                            |
| Volume packaging     | 0.161l                      | 0.36l                       | 1.26l                        | 7.46l                        |
| Dimensions (l/w/h)   | 44 mm/44 mm/106 mm          | 56 mm/56 mm/146 mm          | 90 mm/90 mm/198 mm           | 223 mm/223 mm/191 mm         |
| Layer                | 672                         | 648                         | 192                          | 18                           |
| Pallet               | 2688                        | 1944                        | 576                          | 90                           |
| Category of products | 1                           | 1                           | 1                            | 1                            |
| Customs tariff       | 23099041                    | 23099041                    | 23099041                     | 23099041                     |
| Country of origin    | DE                          | DE                          | DE                           | DE                           |
| Type of packaging    | can                         | can                         | can                          | Eimer                        |

| PU 1 data        |                             |                             |                              |                              |
|------------------|-----------------------------|-----------------------------|------------------------------|------------------------------|
| Product name     | JBL PRONOVO BOTIA TAB 100ml | JBL PRONOVO BOTIA TAB 250ml | JBL PRONOVO BOTIA TAB 1000ml | JBL PRONOVO BOTIA TAB M 5.5l |
| Art. No.         | 3115273                     | 3115300                     | 3115500                      | 3115600                      |
| PU 1 material    | film gr                     | film gr                     | film gr                      | None                         |
| PU 1 weight      | 2.75 g                      | 5.1 g                       | 9.39 g                       | -                            |
| PU 1 length      | 132 mm                      | 168 mm                      | 270 mm                       | -                            |
| PU 1 width/depth | 88 mm                       | 112 mm                      | 90 mm                        | -                            |
| PU 1 height      | 106 mm                      | 146 mm                      | 198 mm                       | -                            |

| PU 2 data        |                             |                             |                              |                              |
|------------------|-----------------------------|-----------------------------|------------------------------|------------------------------|
| Product name     | JBL PRONOVO BOTIA TAB 100ml | JBL PRONOVO BOTIA TAB 250ml | JBL PRONOVO BOTIA TAB 1000ml | JBL PRONOVO BOTIA TAB M 5.5l |
| Art. No.         | 3115273                     | 3115300                     | 3115500                      | 3115600                      |
| PU 2 material    | cardboard gr                | cardboard gr                | cardboard gr                 | None                         |
| PU 2 weight      | 352 g                       | 426 g                       | 468 g                        | -                            |
| PU 2 length      | 284 mm                      | 244 mm                      | 284 mm                       | -                            |
| PU 2 width/depth | 374 mm                      | 354 mm                      | 374 mm                       | -                            |
| PU 2 height      | 223 mm                      | 448 mm                      | 408 mm                       | -                            |

| Trade data        |                             |                             |                              |                              |
|-------------------|-----------------------------|-----------------------------|------------------------------|------------------------------|
| Product name      | JBL PRONOVO BOTIA TAB 100ml | JBL PRONOVO BOTIA TAB 250ml | JBL PRONOVO BOTIA TAB 1000ml | JBL PRONOVO BOTIA TAB M 5.5l |
| Art. No.          | 3115273                     | 3115300                     | 3115500                      | 3115600                      |
| Till receipt text | JBL BOTIA TAB M             | JBL BOTIA TAB M             | BOTIA TAB M                  | JBL BOTIA TAB M              |
| Shelf placement   | -                           | -                           | -                            | -                            |

