





Food tablets for herbivorous sucker catfish (1-20 cm)

Suitable for:









- Food wafer for all sucker catfish that feed on plant and algae, e.g. long-whiskered cat fish. The wood content covers the needs of catfish species which need to grate wood
- Use: allow to sink to the aquarium ground. For twilight and nocturnal species, feed only after switching off the lights
- Natural raw ingredients without artificial additives for a natural diet and healthy aquarium fish. Tested on JBL research expeditions in the tropics
- Very high acceptance even with choosy species thanks to high-quality raw ingredients such as spinach, spirulina and shrimps
- Package contents: original sealed, air- and light-tight, recyclable tin with dosing aid inside the lid (for 100 ml, 250 ml, 1000 ml)



Date: 30.12.2023 Produced by:









Product information

The group of sucker catfish is quite diverse in its diet: the vast majority of species feed herbivorously on plants and algae. Of these, many species are also wood eaters. However, there are also species, such as the zebra or cactus catfish, which do not eat plants but eat micro-organisms, insects and even dead fish. For these species JBL PRONOVO TAB is the right choice. For all other sucking catfish, JBL PRONOVO PLECO has been developed and is available in two sizes.

Some catfish species will not come out when the light is on. If you then feed food tablets, they will disintegrate and tend to be eaten by the other fish, if they like "green food tablets". In this case it is advisable to feed JBL PRONOVO PLECO AFTER switching off the light. The wafers (chips which are thinner than food tablets) do not dissolve so quickly and retain their consistency until the catfish are active, so that the sucker catfish have to grate vigorously (important for the teeth - no joke!).

What do you do if you keep different bottom dwellers in an aquarium? Fish also have a sense of taste and refuse food they don't like. We do not keep ANY bottom-dwelling species, which would experience problems if they ate from a food tablet intended for another species (e.g. loaches). So feed both types of tablets (e.g. JBL PRONOVO PLECO for herbivorous sucker catfish and JBL PRONOVO TAB for carnivorous or omnivorous catfish) at the same time. It is even beneficial for the fish, as they are offered more variety and yet each group of fish gets its species-appropriate nutrition.

| Further information | |
|-----------------------|----------|
| FAQ | ~ |
| Blog | ~ |
| Press | ~ |
| Laboratory/calculator | × |
| Worth reading | ~ |
| Spare parts | × |
| Video | ~ |
| GarantiePlus | × |
| Instructions | ~ |
| QR code | |











Product details

| Article data | | | | |
|--------------------------|------------------------------------|------------------------------------|-------------------------------------|--------------------------------|
| Product name | JBL PRONOVO PLECO WAFER M 100ml | JBL PRONOVO PLECO WAFER M 250ml | JBL PRONOVO PLECO WAFER M 1000ml | JBL PRONOVO PLECO WAFER M 5.5I |
| Art. No. | 3133100 | 3133300 | 3133500 | 3133600 |
| EAN number | 4014162313317 | 4014162313331 | 4014162313355 | 4014162313362 |
| EAN as barcode | | | | |
| Content | 100 ml | 250 ml | 1000 ml | 5,5 |
| Weight | 53 g | 133 g | 530 g | 2900 g |
| Expiry months | 36 | 36 | 36 | 36 |
| RRP incl. VAT | 6,11 € | 11,50 € | 31,38 € | 98,79 € |
| Base price | 115.28 € | 86.47 € | 59.21 € | 34.07 € |
| Nominal filling quantity | 53 g | 133 g | 530 g | 2900 g |
| Base quantity | 1 kg | 1 kg | 1 kg | 1 kg |
| Gross weight | 73.94 g | 168.9 g | 609.96 g | 3079 g |
| Net weight | 53 g | 133 g | 530 g | 2900 g |
| Weight change | - | - | - | - |

| Disposal | | | | |
|-------------------------|------------------------------------|------------------------------------|-------------------------------------|--------------------------------|
| Product name | JBL PRONOVO PLECO WAFER M 100ml | JBL PRONOVO PLECO WAFER M 250ml | JBL PRONOVO PLECO WAFER M 1000ml | JBL PRONOVO PLECO WAFER M 5.5I |
| Art. No. | 3133100 | 3133300 | 3133500 | 3133600 |
| Green dot | ✓ | ✓ | ✓ | ✓ |
| Group electronic waste | - | - | - | - |
| Disposal weight | - | - | - | - |
| Battery type | - | - | - | - |
| Battery return | - | - | - | - |
| Battery rechargeable | - | - | - | - |
| Disposal weight battery | - | - | - | - |
| Non-returnable glass | - | - | - | - |
| PPK | - | - | - | - |
| Plastic small | 20.94 g | 35.9 g | 79.96 g | - |
| Plastic large | - | - | - | 179 g |
| Disposal weight metal | 0 g | 0 g | 0 g | 0 g |

| Features | | | | |
|-------------------|------------------------------------|------------------------------------|-------------------------------------|--------------------------------|
| Product name | JBL PRONOVO PLECO WAFER M 100ml | JBL PRONOVO PLECO WAFER M 250ml | JBL PRONOVO PLECO WAFER M 1000ml | JBL PRONOVO PLECO WAFER M 5.5I |
| Art. No. | 3133100 | 3133300 | 3133500 | 3133600 |
| Animal species | Catfish | Catfish | Catfish | Catfish |
| Animal size | 1 - 20 cm | 1 - 20 cm | 1 - 20 cm | 1 - 20 cm |
| Animal age group | All aquarium fish | All aquarium fish | All aquarium fish | All aquarium fish |
| Volume of habitat | For all aquariums | For all aquariums | For all aquariums | For all aquariums |

Date: 30.12.2023 Produced by:

Country: Brunei Darussalam 3 Language: English









| Features | | | | |
|----------------------|---|---|---|---|
| Material | You will find details of the materials used under Details & Composition | You will find details of the materials used under Details & Composition | You will find details of the materials used under Details & Composition | You will find details of the materials used under Details & Composition |
| Food type | food tablets | food tablets | food tablets | food tablets |
| Colour | green | green | green | green |
| Dosage | Feed 1 – 2 times a day, as much as can be eaten within a few minutes | Feed 1 – 2 times a day, as much as can be eaten within a few minutes | Feed 1 – 2 times a day, as much as can be eaten within a few minutes | Feed 1 – 2 times a day, as much as can be eaten within a few minutes |
| Transport conditions | - | - | - | - |

Date: 30.12.2023 Produced by:









| Electronic label / illuminant | | | | |
|-------------------------------|------------------------------------|------------------------------------|-------------------------------------|--------------------------------|
| Product name | JBL PRONOVO PLECO WAFER M 100ml | JBL PRONOVO PLECO WAFER M 250ml | JBL PRONOVO PLECO WAFER M 1000ml | JBL PRONOVO PLECO WAFER M 5.5I |
| Art. No. | 3133100 | 3133300 | 3133500 | 3133600 |
| Ambient temperature | - | - | - | - |
| Start time | - | - | - | - |
| Mercury | - | - | - | - |
| Tube length | - | - | - | - |
| Service life | - | - | - | - |
| Lumen | - | - | - | - |
| CRI value | - | - | - | - |
| Dimmable | - | - | - | - |
| Switching cycles | - | - | - | - |
| PAR value | - | - | - | - |
| Energy efficiency class | - | - | - | - |
| UV-A | - | - | - | - |
| UV-B | - | - | - | - |
| UV-C | - | - | - | - |
| Colour temperature | - | - | - | - |
| Base designation | - | - | - | - |

| Technical data | | | | |
|-------------------------------|------------------------------------|------------------------------------|-------------------------------------|--------------------------------|
| Product name | JBL PRONOVO PLECO WAFER M 100ml | JBL PRONOVO PLECO WAFER M 250ml | JBL PRONOVO PLECO WAFER M 1000ml | JBL PRONOVO PLECO WAFER M 5.5I |
| Art. No. | 3133100 | 3133300 | 3133500 | 3133600 |
| Range in litres | - | - | - | - |
| Range from - to | - | - | - | - |
| Range in days | - | - | - | - |
| Range tank length | - | - | - | - |
| Output in watts | - | - | - | - |
| Output per hour | - | - | - | - |
| Output per day | - | - | - | - |
| Height | 106 mm | 146 mm | 198 mm | 191 mm |
| Length | 44 mm | 56 mm | 90 mm | 223 mm |
| Width | 44 mm | 56 mm | 90 mm | 223 mm |
| Diameter | 44 mm | 56 mm | 90 mm | 223 mm |
| Voltage | - | - | - | - |
| For | - | - | - | - |
| T8 26mm (watt) | - | - | - | - |
| T5 16mm (watt) | - | - | - | - |
| Size | - | - | - | - |
| Content for | - | - | - | - |
| Filter container volume | - | - | - | - |
| Volume filter media | - | - | - | - |
| Hose connections pressure/out | - | - | - | - |
| Hose connections suction/in | - | - | - | - |

Date: 30.12.2023 Produced by:









| Technical data | | | | |
|----------------|---|---|---|---|
| Delivery head | - | - | - | - |







| Food type | food tablets |
|------------------|--|
| Sub product type | complete food for ornamental fish |
| Dosing | Feed 1 – 2 times a day, as much as can be eaten within a few minutes |

Analytical components



| Crude protein | 29 % |
|---------------|------|
| Fat content | 7 % |
| Crude fibre | 11 % |
| Crude ash | 8 % |

Composition



| wheat germs |
|---------------|
| spinach |
| soy bean meal |
| Wood fibre |
| salmon meal |
| green meal |
| spirulina |
| shrimp meal |
| wheat meal |
| rice meal |
| yeast extract |
| Beet pulp |

Country: Brunei Darussalam

Language: English

il+ Additives

| Colc | urir | σς |
|------|------|----|

are not added by JBL

| Anti | | | |
|------|--|--|--|
| | | | |
| | | | |

1b306(i)

30.12.2023 Date: Produced by:









Additional information for the specialist trade sector

| Article data | | | | |
|----------------------|------------------------------------|------------------------------------|-------------------------------------|--------------------------------|
| Product name | JBL PRONOVO PLECO WAFER M 100ml | JBL PRONOVO PLECO WAFER M 250ml | JBL PRONOVO PLECO WAFER M 1000ml | JBL PRONOVO PLECO WAFER M 5.5I |
| Art. No. | 3133100 | 3133300 | 3133500 | 3133600 |
| VAT | 7% | 7% | 7% | 7% |
| Sales unit (SU) | 6 | 6 | 3 | 1 |
| Volume packaging | 0.1611 | 0.36l | 1.26l | 7.461 |
| Dimensions (l/w/h) | 44 mm/44 mm/106 mm | 56 mm/56 mm/146 mm | 90 mm/90 mm/198 mm | 223 mm/223 mm/191 mm |
| Layer | 672 | 648 | 192 | 18 |
| Pallet | 2688 | 1944 | 576 | 90 |
| Category of products | 1 | 1 | 1 | 1 |
| Customs tariff | 23099041 | 23099041 | 23099041 | 23099041 |
| Country of origin | DE | DE | DE | DE |
| Type of packaging | can | can | can | Eimer |

| PU 1 data | | | | |
|------------------|------------------------------------|------------------------------------|-------------------------------------|--------------------------------|
| Product name | JBL PRONOVO PLECO WAFER M 100ml | JBL PRONOVO PLECO WAFER M 250ml | JBL PRONOVO PLECO WAFER M 1000ml | JBL PRONOVO PLECO WAFER M 5.5I |
| Art. No. | 3133100 | 3133300 | 3133500 | 3133600 |
| PU 1 material | film gr | film gr | film gr | None |
| PU 1 weight | 2.75 g | 5.1 g | 9.39 g | - |
| PU 1 lengh | 132 mm | 168 mm | 270 mm | - |
| PU 1 width/depth | 88 mm | 112 mm | 90 mm | - |
| PU 1 height | 106 mm | 146 mm | 198 mm | - |

| PU 2 data | | | | |
|------------------|------------------------------------|------------------------------------|-------------------------------------|--------------------------------|
| Product name | JBL PRONOVO PLECO WAFER M 100ml | JBL PRONOVO PLECO WAFER M 250ml | JBL PRONOVO PLECO WAFER M 1000ml | JBL PRONOVO PLECO WAFER M 5.5I |
| Art. No. | 3133100 | 3133300 | 3133500 | 3133600 |
| PU 2 material | cardboard gr | cardboard gr | cardboard gr | None |
| PU 2 weight | 352 g | 426 g | 468 g | - |
| PU 2 lengh | 284 mm | 244 mm | 284 mm | - |
| PU 2 width/depth | 374 mm | 354 mm | 374 mm | - |
| PU 2 height | 223 mm | 448 mm | 408 mm | - |

| Trade data | | | | |
|-------------------|------------------------------------|------------------------------------|-------------------------------------|--------------------------------|
| Product name | JBL PRONOVO PLECO WAFER M 100ml | JBL PRONOVO PLECO WAFER M 250ml | JBL PRONOVO PLECO WAFER M 1000ml | JBL PRONOVO PLECO WAFER M 5.5I |
| Art. No. | 3133100 | 3133300 | 3133500 | 3133600 |
| Till receipt text | JBL PLECO WAFER | JBL PLECO WAFER | JBL PLECO WAFER | JBL PLECO WAFER |
| Shelf placement | - | - | - | - |

Date: 30.12.2023 Produced by:

Country: Brunei Darussalam Language: English

